

Navigating Differently: Alternative Drive Controls: Navigating Differently

Jacelyn Goh
Linds Rehabilitation Equipment
Clinical Educator / OT / ATP



The Wheelchair People

Learning Objectives

- Participants will be able to list methods to improve driving skills through standard set up prior to considering alternative driving systems
- Participants will be able to differentiate different driving systems for powered wheelchairs and list the benefits for different client groups
- Participants will be to list 2 considerations for additional functions of alternative driving systems

Alternate Drive Controls – What and Why?

Standard Drive System

- R-Net and Q-Logic
- Proportional control
- With standard joystick / modified joystick
- Mounted on armrest/ other suitable locations
- Expandable vs Non-Expandable
- Function controls on the controller / Easy access



Physical requirements:

- Sufficient UL movement and fine motor control
- Ability to sustain force
- Joystick placement


Q-LOGIC 3
Advanced Drive Control System



Alternate Drive Controls – What and Why?

Adjustability within this standard system:

- Programming (Speed, Acceleration, Active Throw, Dead Band, Tremor dampening, Active Orientation)
- Joystick Handles
- Mounting systems
- Access buttons




Alternate Drive Controls – What and Why?

Alternate Drive Systems

- Different sources of control
 - Placement of systems
 - Would require additional control interface
-
- Why use an alternate drive system?
 - To accommodate for client's needs/ability to drive
 - Better posture / access to support endurance



<https://www.numotion.com/products-services/innovative-lifestyle-products/ability-drive#overview>



The Wheelchair People

Assessment Process for Alternate Drive

Controls:

- Physical Abilities
 - Available Movements (ULs / Head/ Oral)
 - Vision
 - Postural Stability and Displacement
 - Prognosis and impact on abilities
- Wheelchair requirements
 - Power Functions
 - Drive wheel location
 - Communication devices
- Access
 - Environment use
- Cognition / Perception
 - Ability to learn
 - Good safety awareness
 - Problem-solving
 - Spatial perception/awareness
- Carer Requirement
 - Ability to manage complex systems
 - Transfers / daily activities in wheelchair



Alternate Control Drive System - Based on Consistent Movement Patterns

Good / Fair Gross Motor Control, Limited Fine Motor

- **Switch (4) / Proximity Switch**
- Single Switch scanning (single body site)
- CoMove It



JSC for foot drive (FCMH)




JSC for arm drive (ACMH)

Alternate Control Drive System - Based on Consistent Movement Patterns

Good / Fair Gross Motor Control, Limited Fine Motor

- Switch (4) / Proximity Switch
- **Single Switch scanning (single body site)**
- CoMove It




<https://hub.permobil.com/blog/alternative-drive-controls-single-switch-scanning-eye-gaze>

Alternate Control Drive System - Based on Consistent Movement Patterns

Good / Fair Gross Motor Control, Limited Fine Motor

- Switch (4) / Proximity Switch
- Single Switch scanning (single body site)
- **CoMove It**



Alternate Control Drive System - Based on Consistent Movement Patterns

Good Fine motor, Limited Gross motor

- Mini and Micro- Switches
 - Significantly lesser amount of force for displacement (~10-50g)
 - Multiple mounting options



Alternate Control Drive System - Based on Consistent Movement Patterns

No UL, Good Head Control

- **Chin control (micro switches)**
- Gyro systems (Vigo and Munevo Glasses)
- Head Controls (proportional vs Head array)



Alternate Control Drive System - Based on Consistent Movement Patterns

No UL, Good Head Control

- Chin control (micro switches)
- **Gyro systems (Vigo and Munevo Glasses)**
- Head Controls (proportional vs Head array)



<https://www.sunrisemedical.com/power-wheelchairs/switch-it-electronics/head-controls/vigo>



The Wheelchair People

Alternate Control Drive System - Based on Consistent Movement Patterns

No UL, Good Head Control

- Chin control (micro switches)
- Gyro systems (Vigo and Munevo Glasses)
- **Head Controls - Switched Head array (Permobil total control, Switch-it Duo Pro, ASL Fusion, Stealth i-Drive)**



 **Dual Pro™**



SWITCH-IT
SMART TECHNOLOGIES



Stealth Products
i-Drive® Head Arrays



The Wheelchair People

Alternate Control Drive System - Based on Consistent Movement Patterns



Fair Head Control, Limited Extremity control

- **Switched Head Array**
- Munevo Glasses

Considerations:

- Head supports required
- Endurance
- Increased adjustments to profiles to reduce postural displacement

Alternate Control Drive System - Based on Consistent Movement Patterns

**Little head control,
Good oral control**

- Sip n Puff,
Combination



Alternate Control Drive System - Based on Consistent Movement Patterns

**Minimum head
control, good eye gaze
control**
- Eye Gaze



The Wheelchair People

Mounting Systems

Setting up the most ergonomic position of drive method to ensure endurance and ease of use

Considerations for carer staff:

- Ease to remove/replace
- Adjustments required
- Durability

Mounting Systems



Tray mounted Systems



Arm rest
mounts



Center Drive System / Mounts



The Wheelchair People

Chin Control / Head Array mounts



Recommendation

“The factors identified in research related to non-use of provided assistive technology (as listed) should be considered by the therapist during the wheelchair prescription process, as these may influence the outcomes.”

https://www.aci.health.nsw.gov.au/__data/assets/pdf_file/0003/167286/Guidelines-on-Wheelchair-Prescription.pdf Page 25





Questions?
Please come and try out the
systems!

Jace Goh
Occupational Therapist/ Product Consultant
jacelyn@lindsrehab.com.au

The Wheelchair People