

Over Toilet Frame / Commode



This versatile product can be used as a freestanding, over toilet frame or a standalone commode. With a maximum user capacity of 130 kg, the height adjustable steel frame features padded armrests for added comfort and anti-slip rubber tips for stable grip, even on wet surfaces. A simple and secure push button mechanism allows the height of each leg to be adjusted to suit individual user requirements.

The plastic contoured seat includes pan and lid, which slides out to facilitate easy and hygienic cleaning. Designed to fit the majority of toilets on the market.

CODE	ED1249
RAP	BE06
WEIGHT	6.6 kg
OVERALL HEIGHT	605 > 760 mm
SEAT HEIGHT	405 > 560 mm
ARMREST HEIGHT	200 mm above seat
OVERALL WIDTH	530 mm
BETWEEN ARMRESTS	460 mm
CONFORMS TO	ARTG 176831

