



VOGUE WALKER RANGE

aspire for...

Comfort
Stability
Style



www.aspirecare.com.au

1300 133 120



VOGUE CARBON FIBRE SEAT WALKER



aspire
vogue

One of the lightest seat walkers on the market! Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super lightweight carbon fibre frame for a superior walking experience.

Signature
Carbon Fibre
Weave



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



Carbon - Black



Carbon -
Candy Red



AVAILABLE
IN 2x SIZES:
MEDIUM
& TALL



Carbon Fibre Frame

Extremely lightweight to
maximise portability

LESS THAN
5KG*

MODEL

WAF705350 | Vogue Carbon Fibre - Medium - Black
WAF705360 | Vogue Carbon Fibre - Tall - Candy Red

FEATURES & BENEFITS shared by the Vogue Carbon Fibre and Lightweight Aluminium Seat Walkers



European Design

Ergonomic design offers
personalised height adjustment,
promoting correct posture



Compact

Turning circle for
improved mobility in tight
spaces



Simple

Folding mechanism locks
for easy transportation and
storage



Handbrakes

Lock for user safety
and stability

	LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT** (kg)	BAG S.W.L (kg)
MED	680	600	695 - 920	240	460	550	225	8"	4.9	6.2	5
TALL	680	605	770 - 1040	240	460	620	225	8"	5.2	6.3	5

VOGUE LIGHTWEIGHT SEAT WALKER



aspire
vogue

This modern seat walker has a clean and crisp frame design, making it a standout when going for a walk or browsing the shops!



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



Red



Champagne

WEIGHS
JUST OVER
5KG



Discreet Cable
Management



**Lightweight
Aluminium Body**

Improves manoeuvrability
and portability

MODEL

WAF705310RE | Vogue Lightweight Aluminium - Red
WAF705310CH | Vogue Lightweight Aluminium - Champagne



Storage

Removable bag
carries up to 5kg
of essentials



Backrest

Deep seat and wide
backrest for improved
comfort and support



Bonus

Walking cane holder

**1 YEAR LIMITED
WARRANTY**

✓ AS/NZ Approved

✓ Vogue is registered in the
ARTG, AUST R ARTG 271255

OVERALL			SEAT			FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT** (kg)	BAG S.W.L (kg)
LENGTH (mm)	WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	WIDTH (mm)	HEIGHT (mm)					
655	620	785 - 920	240	465	540	220	8"	5.3	6.4	5

* Carbon Fibre Walker - Medium

** Includes accessories, bag and backstrap.

VOGUE ADVENTURE SEAT WALKER



aspire
vogue

Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super sized front wheels can manage any surface, perfect for the adventurous user!

Wider Seat
supports up to
200kg



Folds to
245mm



Large
Multi-Terrain
Wheels



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



Black

AutoLocate

Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate.



MODEL

WAF705450 | Vogue Adventure - Black

FEATURES & BENEFITS



European Design

Ergonomic design offer personalised height adjustment, promoting correct posture



Compact

Turning circle for improved mobility in tight spaces



Simple

Folding mechanism locks for easy transportation and storage

1 YEAR LIMITED WARRANTY

✓ AS/NZ Approved

✓ Vogue is registered in the ARTG, AUST R ARTG 271255



Handbrakes

Lock for user safety and stability



Storage

Removable bag carries up to 5kg of essentials



Bonus

Walking cane holder

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)
770	695	900 - 1035	230	540	530	245	Front 11" Rear 9"	8.8	10	5

VOGUE FOREARM WALKER



For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



White / Black

**Lightweight
Aluminium Body**
Improves portability
and manoeuvrability



MODEL

WAF705400 | Vogue Forearm - White/Black

FEATURES & BENEFITS



European Design
Ergonomic design offer personalised height adjustment, promoting correct posture




Simple
Folding mechanism locks for easy transportation and storage



Bonus
Walking cane holder



Handbrakes
Lock for user safety and stability



Storage
Removable bag carries up to 5kg of essentials

1 YEAR LIMITED WARRANTY

✓ AS/NZ Approved

✓ Vogue is registered in the ARTG, AUST R ARTG 271255

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)
660	620	1020 - 1270	240	460	620	225	8"	9.1	10*	5

* Includes accessories, bag and backstrap.

VOGUE INDOOR WALKER



aspire
vogue

Ideal for use within the home, the Indoor Walker is perfect for portable dining and in home mobility assistance.



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



Red / Black



Champagne
/ Black



Folds to
270mm



Portable Dining

Safely collect and transport
meals or hot drinks

On-cart Dining



Turn Around.
Slide Over Knee.
Lock & Enjoy.



MODEL

WAF705100RE | Vogue Indoor - Red/Black
WAF705100CH | Vogue Indoor - Champagne/Black


FEATURES & BENEFITS



European Design
Ergonomic design offer
personalised height
adjustment, promoting
correct posture



Compact
Turning circle for
improved mobility in
tight spaces



Conveniently
Stands up when
folded for flat
storage



Handbrakes
Lock for user safety
and stability



Storage
Removable bag carries
up to 5kg of essentials



Bonus
Walking cane
holder

1 YEAR LIMITED WARRANTY

✓ AS/NZ Approved

✓ Vogue is registered in the
ARTG, AUST R ARTG 271255

OVERALL			FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)	TRAY S.W.L (kg)
LENGTH (mm)	WIDTH (mm)	HEIGHT (mm)						
650	540	860 - 970	270	6"	4.4	5.1	5	5

* Includes accessories, bag and tray.



Oxygen
Bottle Holder Bag
WAA695600



Night Light - LED Lamp
WAA695630



Single Hand Brake System* -
Compatible with Aspire Vogue Carbon Walkers
WAA695640

**Single Hand Brake System are not compatible with the Vogue Indoor Walker*



CRUTCHES, FRAMES & STICKS



SCOOTERS



OVERBED TABLES



BED STICKS & POLES



RAMPS



AIR MATTRESSES & FOAM MATTRESSES



COMMODOES, URINALS & BED PANS



OVER TOILET AIDS



SHOWER STOOLS & CHAIRS



RAISED TOILET SEATS



SHOWER COMMODOES



BATH BOARDS & BENCHES



ADJUSTABLE BEDS



POWER WHEELCHAIRS



PRESSURE CUSHIONS



LIFT & RECLINE CHAIRS



Distributed by:

power mobility
The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034
sales@powermobility.com.au • powermobility.com.au

ASBR21.43 130820

