

J 623 2.0

Package Inclusions

- Quick Height Adjustable (QHA) armrests
- TB3 Tilt - TRU-Balance 50 degree tilt
- Swingaway legrests with angle adjustable foot plates or centre-mount footboard
- Choice of either TRU-Comfort back or upgrade to ADI backs
- Choice of cushion
- Stealth headrest

Further options available - refer to script form.

With features including powerful motors that achieve speeds up to 10 km/h and improved responsiveness, the Jazzy® 623 2.0 is engineered to meet the performance needs of the most active user. The Jazzy® 623 2.0 comes standard with TRU-Balance 3 50 degree tilt and optional 25.4 cm (10") of lift elevation, while a VR2 controller offers two seating functions (as standard) and delivers smooth, nimble manoeuvrability matched with exceptional outdoor performance.

Standard Colours



Alternative colour options available - refer to script form.



Attendant Control



Jazzy 623 2.0 shown in red with activated tilt



QHA Armrests



Lift



Specifications

DRIVE WHEELS	Standard: 35.56 cm (14") Solid Optional: 35.56 cm (14") Pneumatic
CASTOR WHEELS	Front: 15.24 cm (6") Solid (standard) 20 cm (8") Solid (optional) Rear: 15.24 cm (6") Solid (standard) 20 cm (8") Solid (optional)
	N/A
MAXIMUM SPEED ⁴	Up to 9.6 km/h (6 mph)
GROUND CLEARANCE	6.9 cm (2.7") (battery tray)
TURNING RADIUS ⁵	52.1 cm (20.5")
OVERALL LENGTH	90.2 cm (35.5") without foot riggings
BASE WIDTH	62.2 cm (24.4")
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 30 cm - 56 cm (12-22") Depth: 30 cm - 56 cm (12-22")
TRU-Balance® 3 Power Tilt	Width: 30.5 cm - 56 cm (12-22") Depth: 30.5 cm - 56 cm (12-22")
SEAT-TO-FLOOR HEIGHTS	
Synergy®/TRU-Balance® 3 Static Seat ⁸	42 cm - 47 cm (16.5-18.5")
TRU-Balance® 3 Power Tilt ⁸	44 cm - 49.5 cm (17.5-19.5")
TRU-Balance® 3 Lift & Tilt ⁸	44 cm - 49.5 cm (17.5-19.5")
MANUAL TILT/RECLINE ¹⁰	Yes
BATTERY SIZE ⁶	NF-22, 55 AH
BATTERY WEIGHT ¹⁴	17.24 kg (38 lbs.)
AVAILABLE ELECTRONICS ¹²	VR2 Controller
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole Mid-Wheel 6®
WEIGHT CAPACITY ⁷	136 kg (300 lbs.) ⁹
BASE WEIGHT	63 kg (138.6 lbs.)



Swingaway legrests with angle adjustable foot plates or centre-mount footboard



Stealth headrest

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems
- 7) Weight capacity is dictated by the selection of seating
- 8) All measurements +/- .3175 cm (1/8") with 46 cm x 46 cm (18" x 18") seat, fixed towers, 136 kg (300 lbs.) weight capacity and solid tyres (add .95 cm (3/8") with pneumatic tyres). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified
- 9) 113 kg (250 lbs.) weight capacity with TRU-Balance® Systems
- 10) The Synergy® Manual Tilt weight capacity is 34 kg (75 lbs.) and maximum dimensions are 35.5 cm x 35.5 cm (14" W x 14" D). Maximum dimensions of the Manual Recline vary by power base and seating system
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 1.4 kg (3 lbs.) based on manufacturer
- 15) Call for availability

