

# Simply Fold and Travel



**Jazzy®**  
**PASSPORT**



Experience the beauty of **lightweight travel** with the **Jazzy® Passport!** This compact power chair folds in just a few simple steps for effortless transport. With convenient storage, the Jazzy® Passport makes traveling a breeze.

Folded Size  
78.7 x 40.6 x 60 cm  
(31" x 16" x 23.5")

PRIDE® **MOBILITY SCOOTERS**  
LIVE YOUR BEST®

[www.pridemobility.com.au](http://www.pridemobility.com.au)



## Specifications

Weight Capacity	113 kg (250 lbs.)
Maximum Speed <sup>1</sup>	Up to 5.8 km/h (3.6 mph)
Ground Clearance <sup>3</sup>	7.9 cm (3") at motor
Turning Radius <sup>3</sup>	61.6 cm (24.25")
Overall Length <sup>3</sup>	101.6 cm (40")
Overall Width <sup>3</sup>	59.7 cm (23.5")
Seat-to-Ground Height	48.3 cm (19")
Seat-to-Foot Platform Height	40 cm (15.74") ( <i>no adjustment</i> )
Casters	20.3 cm (8") solid
Drive Tires	30.5 cm (12") Pneumatic drive
Suspension	Limited
Range Per Charge <sup>1,2</sup>	Up to 15 km (9.6 miles)
Total Weight Without Batteries <sup>3, 5</sup>	24 kg (54 lbs.)
Heaviest Piece When Disassembled <sup>3</sup>	N/A
Standard Seating	Material: Gray, Foam Dimensions: 45 x 43 cm (17.7" x 16.9")
Drivetrain	2x 180W brushed motors
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>3</sup>	25.2V 18.2Ah Lithium ion
Battery Weight <sup>4</sup>	4 kg (8.81 lbs.)
Battery Charger	2-amp, off-board
Colour	 Silver

## Features

- A compact and lightweight design makes traveling a breeze
- Back seat pocket for extra storage while on-the-go
- Folds in just a few simple steps for effortless transport
- Includes convenient, under seat mesh storage bag, mesh cup holder and 101.6 cm (40") lap belt
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extreme comfort

## Accessories

- XLR USB charger



- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

**power: mobility**  
The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034

sales@powermobility.com.au • powermobility.com.au



20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)



**Pride<sup>®</sup>**  
Mobility Products  
Australia Pty. Ltd.