

STELLAR GL



The ***Stellar GL*** is a versatile 45° low knee-rise tilt chair. It promotes self propulsion, table access, and ***participation in mobility related activities of daily living***. The rigid frame provides superior stability and maneuverability.



Hand
Propulsion



Foot
Propulsion



Self-Tilt



Desk/Table
Accessible



Transit Tie
Downs

Exclusive features

- ✓ Foot and hand propel in tilt position with low knee-rise
- ✓ Short wheelbase for maneuverability
- ✓ Simple center of gravity adjustments

Available options See order form

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt



Stellar GLT 45° tilt chair with a lightweight aluminum frame

Stellar GL 45° tilt chair with a rugged steel frame

The lightweight design of the ***Stellar GLT*** combines the benefits of 45° tilt with a high degree of mobility independence.

YOUR PICK

The rugged steel frame of the standard ***Stellar GL*** makes it a great choice for heavier users or those with high tone/agitation.

STELLAR GL



Stellar GLT



DYNAMIC TILT
Up to 45°



WEIGHT CAPACITY
250 lbs



SEAT WIDTH
14" to 22"



SEAT DEPTH
16" to 22"



SEAT HEIGHT
13" to 20"



BACK SEAT ANGLE
90° to 122°



OVERALL WIDTH
+9.5" to 10.5"



CHAIR WEIGHT
34 lbs

visit pdgmobility.com/glt
for more information



Stellar GL

visit pdgmobility.com/gl
for more information



DYNAMIC TILT
Up to 45°



WEIGHT CAPACITY
Up to 450 lbs



SEAT WIDTH
14" to 24"



SEAT DEPTH
16" to 22"



SEAT HEIGHT
13" to 20"



BACK SEAT ANGLE
90° to 122°



OVERALL WIDTH
+9.5" to 10.5"



CHAIR WEIGHT
51 lbs

