

Jazzy Air²



Shown in Pearl White

136 kg

Weight Capacity

Up To

6.4 km/h

Maximum Speed

Elevated & Non-Elevated

61.5 cm

Turning Radius

Up To

27 km

Range Per Charge*



Bright LED marker lights

At the touch of a switch, the **Jazzy Air²** elevates 30.4 cm in just 11 seconds. With a sleek, more modern look, this power chair safely drives at speeds up to 6.4 km/h while elevated, allowing users to enjoy face-to-face social engagement. With Mid-Wheel 6[®] Drive and Active-Trac[®] Suspension, the Jazzy Air 2 offers superb performance and stability both indoors and out.

Specifications

| | |
|-------------------------------------|---|
| Model Name | Jazzy Air® 2 |
| Weight Capacity | 136 kg (300 lbs.) |
| Maximum Speed ¹ | 6.4 km/h (4 mph) |
| Ground Clearance ² | 2.54 cm (1") at motor bracket |
| Turning Radius ¹ | 61.5 cm (24.25") |
| Overall Length ² | 106.7 cm (42") |
| Overall Width ² | 64.1 cm (25.25") |
| Front Wheels | 15.24 cm (6") flat free |
| Drive Wheels | 25.4 cm (10") flat free |
| Rear Wheels | 15.24 cm (6") flat free |
| Drivetrain | Two motor, Mid-Wheel 6® |
| Braking System | Regenerative and electromechanical |
| Suspension Type | Active-Trac (ATX) with Mid-Wheel 6® Technology |
| High-Back Seat Size | 45.7 x 45.7 - 50.8 cm (18" x 18" - 20") |
| Seat-to-Floor Range | 54 cm - 83.8 cm (21.25" - 33") |
| Range Per Charge ^{2,4} | Up to 27 km (16.8 miles) |
| Standard Electronics | 40A, PG VR2 |
| Component Weights ⁶ | Base: 40 AH Batteries 120 kg 40 AH Batteries: 14.28 kg |
| Battery Requirements ^{3,5} | (2) 12V, 40 AH |
| Battery Charger | Off-board, 3.5 amp |
| Warranty | Frame: Lifetime limited Drivetrain: 2-year Electronics: 1-year Battery: (Gel battery) 1-year (AGM battery) 6-months |

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. AGM or gel cell type required.
4. Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on battery manufacturer.
6. Weight includes base, seat, and batteries.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition, and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius, and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034
sales@powermobility.com.au
powermobility.com.au

Elevate your life to the next level!

- Newly re-designed lift features 30.4 cm of power adjustable seat height in 11 seconds
- Convenient "Air" switch on armrest activates elevation
- Active-Trac® Suspension with Mid-Wheel 6® Technology allows for maximum stability in the elevated and non-elevated positions
- Bright LED marker lights offers visibility in all environments
- High-back reclining sport seat
- Height and angle adjustable headrest
- Front under-seat storage keeps your valuables accessible

Additional Features

- Angle, height and depth adjustable armrests
- Integrated foot platform with rubber insert and 5 cm of adjustment for comfort in both the elevated and non-elevated positions
- Single-sided castor forks for enhanced manoeuvrability


Available Accessories

- Rear bag / walking stick holder**
- Walker holder***
- Oxygen tank holder***
- Mobile phone holder
- Cup holder
- Saddle Bag**
- Under-seat storage compartment

*Must have Pride rear bag to attach

**Requires supporting bracket

Available Colours



20-24 Apollo Drive, Hallam VIC 3803
ACN 088 609 661 ABN 66 088 609 661

AU_JAZZYAIR2_01-09-20