

# Jazzy®

## ELITE HD

The Jazzy® Elite HD features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 35.6 cm (14") drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.



*Jazzy® Elite HD  
shown in Jazzy Red  
with high-back seat*

## SPECIFICATIONS

<b>Weight Capacity</b>	204 kg (450 lbs.)*
<b>Maximum Speed</b>	Up to 6.4 km/h (4 mph)
<b>Ground Clearance</b>	7.6 cm (3")
<b>Turning Radius</b>	61 cm (24")
<b>Overall Length</b>	90.2 cm (35.5")
<b>Overall Width</b>	61.6 cm (24.25")
<b>Drive Wheels</b>	35.6 cm (14") solid
<b>Front Wheels</b>	7.6 cm (3") solid anti-tips
<b>Rear Wheels</b>	15.2 cm (6") solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, front-wheel drive
<b>Braking System</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard Electronics</b>	60 amp, Dynamic Shark controller
<b>Battery Charger</b>	Off-board, 5 amp
<b>Per-Charge Range</b>	Up to 15.77 km (9.8 miles)
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, NF-22
<b>Battery Weight</b>	17 kg (37 lbs.) each
<b>Base Weight</b>	44 kg (98.5 lbs.)
<b>Standard Seat Weight</b>	23 kg (51.5 lbs.) (high-back)
<b>Standard Seat</b>	51 x 46-51 cm (20" x 18-20")

## FEATURES

- Two-motor, front-wheel drive
- 35.5 cm (14") knobby tyres and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- 5 amp, off-board charger

## OPTIONS & ACCESSORIES

- Cup holder
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag



Elite HD shown in Jazzy Blue

## COLOUR CHOICES



Jazzy Red



Jazzy Blue

## OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

**power: mobility**  
The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034  
sales@powermobility.com.au  
powermobility.com.au

**Pride**  
Mobility Products  
Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)



The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

\* Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.



AU-JZEIteHD\_083017