

Single Fold RAMPS

PRIDE **R** MPS

Easy single-fold design



272kg
Weight Capacity

75cm
Wide

Portable
with a High-Traction Surface



Single Fold Ramps assist in traversing curbs, steps, & more

Single Fold Ramps from Pride Mobility® offer a reliable solution for transitioning curbs, thresholds, and other obstacles. They are available in several different lengths and can accommodate a weight capacity of 272 kg (600 lbs.). Moreover, lightweight Single Fold Ramps are easily portable, making them ideal for home and travel.

For more information visit pridemobility.com.au or call **03 8770 9600**
20-24 Apollo Drive, Hallam VIC 3803 | ACN 088 609 661 | ABN 66 088 609 66


PRIDE
MOBILITY
LIVE YOUR BEST®



Features

- Handle for easy transport
- Lightweight & portable
- Integrated, high-traction surface
- Competitively priced
- 75 cm (29.5") wide when folded
- 272 kg (600 lbs.) Weight capacity
- Available in the following sizes:

Size	Weight
61 cm (2')	5.8 kg (12.78 lbs)
91 cm (3')	7.7 kg (16.97 lbs)
122 cm (4')	10.3 kg (22.7 lbs)
152 cm (5')	2.8 kg (28.22 lbs)
183 cm (6')	15.4 kg (33.95 lbs)



Ramp Length	Rise				
	16.1° (29%)	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)
61 cm (2')	17° (30%)	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)
91 cm (3')	17° (30%)	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)
122 cm (4')	17° (30%)	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)
152 cm (5')	17° (30%)	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)
183 cm (6')	17° (30%)	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)
	16.1° (29%)	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)
	14.5° (26%)	12.8° (23%)	11.2° (20%)	9.6° (17%)	8° (14%)
	12.8° (23%)	11.2° (20%)	9.6° (17%)	8° (14%)	7.2° (13%)
	11.2° (20%)	9.6° (17%)	8° (14%)	7.2° (13%)	6.4° (11%)
	9.6° (17%)	8° (14%)	7.2° (13%)	6.4° (11%)	5.6° (10%)
	8° (14%)	7.2° (13%)	6.4° (11%)	5.6° (10%)	4.8° (8%)
	7.2° (13%)	6.4° (11%)	5.6° (10%)	4.8° (8%)	4° (7%)
	6.4° (11%)	5.6° (10%)	4.8° (8%)	4° (7%)	3.8° (7%)
	5.6° (10%)	4.8° (8%)	4° (7%)	3.8° (7%)	3.2° (6%)
	4.8° (8%)	4° (7%)	3.8° (7%)	3.2° (6%)	2.9° (5%)
	4° (7%)	3.8° (7%)	3.2° (6%)	2.9° (5%)	2.4° (4%)
	3.8° (7%)	3.2° (6%)	2.9° (5%)	2.4° (4%)	1.9° (3%)
	3.2° (6%)	2.9° (5%)	2.4° (4%)	1.9° (3%)	-
	2.9° (5%)	2.4° (4%)	1.9° (3%)	-	-
	2.4° (4%)	1.9° (3%)	-	-	-
	1.9° (3%)	-	-	-	-
	-	-	-	-	-
	-	-	-	-	2.5 cm (1")