

# PRIDE PUID RECLINER COLLECTION

COLLECTION





Pride® Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance and are available in several beautiful fabrics and styles to fit any décor. The state-of-the-art lift system utilised in Pride Lift Chairs will recline you further back than any other lift chair. Plus, Pride Lift Chairs are equipped with an emergency battery backup system fail in your home.

### LIFT CHAIR **DESCRIPTIONS**

**C1** 

3-Position Lift Chair

3-Position Chaise Lounger

Infnite-Position Chaise Lounger

C-101

3-Position Chaise Lounger

**LC-101 Euro Leather** 

3-Position Chaise Lounger (With European Style Rounded Arm Design)

**D30** 

3-Position Chaise Lounger

3-Position Chaise Lounger

Infnite-Position Chaise Lounger

Infnite-Position Chaise Lounger, Featuring Trendelenburg Position (Elevates Legs Above Head for Greater Relaxation) & Zero-Gravity Position

**LC-358XL** 

Bariatric 3-Position Chaise Lounger

LC-107

Infnite-Position Chaise Lounger (With Extra Leg Extension Length)

## AVAILABLE FABRICS 1



Ocean



**Burgundy** 





Latte





Sandal

**LEATHER** 



Red



**LEVANTE** VINYL



**Burgundy** 



**Royal Blue** 

<sup>1.</sup> Not all fabrics are available for every model. Please see Specifications for details.

<sup>2.</sup> This chair does not recline.

# LIFT CHAIR MODELS









Accented Pillow-back Design Shown in Arctic Blue

Sewn Pillow-back Design Shown in Arctic Blue Sewn Pillow-back Design Shown in Chocolate Button-back Design Shown in Sandal









LC-101 Euro Leather

Split-back Design
Shown in Beige Leather

D30

Split-back Design Shown in Burgundy

Pillow-back Design Shown in Burgundy 660
Pillow-back Design
Shown in Latte







Mongolian-back Design Shown in Black (Leather)

Split-back Design Shown in Pecan

Button-back Design Shown in Latte

LC-358XL





PLEASE NOTE: 2-position is a TV recline.
3-position is a full recline.
Infinite-position is a full recline with independent leg rest function.

### **SPECIFICATIONS**

	C1	<b>C5</b> <sup>2</sup>	C6 <sup>2</sup>	C-101	LC-101 Euro Leather	D30
Back Width:	51 cm (20")	71 cm (28")	71 cm (28")	72 cm (28.5")	58 cm (23")	74 cm (29")
Seat-To-Floor:	41 cm (16")	46 cm (18")	48 cm (19")	47 cm (18.5")	47 cm (18.5")	48 cm (19")
Seat Depth:	51 cm (20")	46 - 51 cm (18 - 20")	53 cm (21")	51 cm (20")	51 cm (20")	56 cm (22")
Seat Width:	46 cm (18")	56 cm (22")	56 cm (22")	51 cm (20")	51 cm (20")1	56 cm (22")
Top of Back to Seat:	66 cm (26")	66 - 71 cm (26 - 28")	66 - 71 cm (26-28")	70 cm (27.5")	70 cm (27.5")	66 cm (26")
Weight Capacity:	113 kg (250 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	170 kg (375 lbs.)
Colour Options:	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal <b>Levante Vinyl:</b> Royal Blue	Cocoa Fern Sandal	<b>Leather:</b> Beige Brown Black Red	Ocean Burgundy Latte Levante Vinyl: Burgundy Royal Blue

	560	660	T-3	LC-107	LC-358XL
Back Width:	71 cm (28")	71 cm (28")	69 cm (27")	74 cm (29")	91 cm (36")
Seat-To-Floor:	48 cm (19")	48 cm (19")	48 cm (19")	51 cm (20")	51 cm (20")
Seat Depth:	48 cm (19")	51 cm (20")	51 cm (20")	53 cm (21")	56 cm (22")
Seat Width:	48 cm (19")	48 cm (19")	51 cm (20")	51 cm (20")	66 cm (26")
Top of Back to Seat:	63.5 cm (25")	63.5 cm (25")	68.5 cm (27")	70 cm (27.5")	74 cm (29")
Weight Capacity:	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	147 kg (325 lbs.)	227 kg (500 lbs.)
Colour Options:	Ocean	Latte	Black (Leather)	Pecan	Burgundy (Vinyl)
	Burgundy			Ocean	Latte
				Burgundy	

#### **BENEFITS**

- Engineered furniture laminate/hardwood frame
- Patented, quiet & smooth lift system
- Large side pocket <sup>1</sup>

- Hand control can be easily moved to left-hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup



# 07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034 sales@powermobility.com.au powermobility.com.au

Note: Measurements may vary +/- 2.5 cm (1"). Specifications subject to change without prior notice.

- Large side pocket not available on LC-440.
   C5 and C6 also available in Pride®'s custom fabric selections range. Please see "Custom Fabric Selections" brochure for further details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



20-24 Apollo Drive Hallam, Victoria 3803 www.pridemobility.com.au ACN 088 609 661 ABN 66 088 609 661



Tel: 03 8770 9600

Fax: 03 9703 2720