

PRIDE[®]
MOBILITY



<27 km

Range per Charge



136 kg

Weight Capacity



61.5 cm

Turning Radius



Jazzy
air₂[®]

Experience Living Your Life at New Heights

Elevate your life to the next level with the Jazzy[®]Air 2. Featuring a captain's seat that elevates up to 30.4 cm, it allows you to rise to eye level during conversations and move at up to 6.4 km/h in the elevated position, keeping you shoulder to shoulder with your peers throughout your day. Everyday tasks become more accessible allowing you to effortlessly extend your reach to counter tops, cabinets, and upper shelves. With Mid-Wheel 6[®] Drive and Active-Trac[®] Suspension, the Jazzy[®]Air 2 offers enhanced performance and stability both indoors and out.

Visit pridemobility.com.au for more information



Jazzy[®] air₂



Narrow Base Width of 55.8 cm



Under Seat Storage



Mid-Wheel 6[®] Technology

RANGE PER CHARGE (MAX) ¹	Up to 27 km
MAXIMUM WEIGHT CAPACITY	136 kg
TURNING RADIUS	61.5 cm
MAXIMUM SPEED ¹	6.4 km/h
LENGTH & WIDTH ²	L: 106.7 cm W: 64 cm
WEIGHT OF HEAVIEST PIECE ⁴	91.6 kg
GROUND CLEARANCE	3.6 cm at motors
SEAT SIZE	W: 40.6 cm x D: 45.7 cm (16" x 18" - 20") W: 45.7 cm x D: 45.7 cm (18" x 18" - 20") W: 50.8 cm x D: 45.7 cm (20" x 18" - 20")
18" X 18" STANDARD SEAT WEIGHT	23 kg (high-back)
HIGH-BACK SEAT-TO-FLOOR RANGE	54 cm - 83 cm

STANDARD ELECTRONICS	40A, PG VR2
DRIVE TRAIN	Two motor, Mid-Wheel 6 [®]
BRAKING SYSTEM	Regenerative and electro-mechanical
SUSPENSION TYPE	Active-Trac (ATX)
FRONT CASTORS	6" solid
DRIVE WHEELS	10" solid
REAR CASTORS	6" solid
BATTERY REQUIREMENTS ³	(2) 40Ah, 12V
BATTERY CHARGER	Off-board, 3.5A
BATTERY WEIGHT	14.3 kg each
COLOUR OPTIONS	<input checked="" type="checkbox"/> Black <input type="checkbox"/> Pearl White <input checked="" type="checkbox"/> Garnet Red <input checked="" type="checkbox"/> Pink Topaz

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Deep cycle (AGM or gel cell type recommended)
4. Seat weight included, Battery weight not included.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

power mobility
The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034
sales@powermobility.com.au • powermobility.com.au

PRIDE[®]
MOBILITY

20-24 Apollo Drive, Hallam, Victoria 38031
(03) 8770 9600 | www.pridemobility.com.au

