

Q6 EDGE 3 *Stretto*™



Q6 EDGE 3 STRETTO™
shown with iLevel® power seat elevation,
Q-Logic 3 environmental controls, and
Stealth® TRU-Comfort 2 cushion.



The Q6 Edge 3 Stretto™ features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 31.8 cm or 35.6 cm drive wheels and optional 30.5 cm of iLevel® power adjustable seat height.



Narrowest
complex
rehab base
on the market


QUANTUM®
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	35.6 cm (14") or 31.8 cm (12.5") pneumatic or flat free tyres
CASTER WHEELS	Front: 12.7 cm (5") flat free Rear: 15.2 cm (5.9") flat free
MAXIMUM SPEED ¹	Up to 10 km/h Up to 5.6 km/h w/iLevel activated
GROUND CLEARANCE	7.4 cm (2.9")
TURNING RADIUS ²	53.3 cm (20.9")
OVERALL LENGTH ²	82.5 cm (32.4") without front riggings
BASE WIDTH	55.2 cm (21.7")/35.6 cm (14") d/wheels 52.7 cm (20.7")/31.8 cm (12.5") d/wheels
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 30.5 cm (12") - 51 cm (20") Depth: 30.5 cm (12") - 51 cm (20")
TRU-Balance® 3 Power Tilt	Width: 30.5 cm (12") - 51 cm (20") Depth: 30.5 cm (12") - 51 cm (20")
TRU-Balance® 3 Tilt and iLevel®	Width: 30.5 cm (12") - 51 cm (20") Depth: 30.5 cm (12") - 51 cm (20")
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	40.6 cm (15.9") - 45.7 cm (17.9")
TRU-Balance® 3 Power Tilt ⁵	44.5 cm (17.5") - 49.5 cm (19.4")
TRU-Balance® 3 Tilt & iLevel® ⁵	44.5 cm (17.5") - 49.5 cm (19.4")
MANUAL RECLINE ⁶	Yes (TRU-Balance® Flex only)
BATTERY SIZE ³	55 AH (2)
BATTERY WEIGHT ⁷	55 AH: 18.1 kg each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6®
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT WITHOUT SEATING SYSTEM	57.3 kg

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

2) Without front riggings

3) Battery capacity may change when adding power positioning systems.

4) Weight capacity is dictated by the selection of seating.

5) All measurements +/- 0.3 cm with 45.7 cm x 45.7 cm seat, fixed towers, 136 kg weight capacity and solid tires (add 0.95 cm with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.

7) Battery weight may vary +/- 1.4 kg based on manufacturer.

8) All Quantum Rehab® units can be ordered with specialty controls.

9) Call for availability.


Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 30.5 cm (12") of seat elevation at walking speed (up to 5.6 km/h)
- Fender lights and USB charger included


Options

- iLevel® option allows 30.5 cm (12") of seat elevation while at walking speed (5.6 km/h)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- Trilogy vent tray⁹
- Compact vent tray⁹
- 30.5 cm (12") x 30.5 cm (12") seat sizes with iLevel (SOF only)

Colours*



RED
CORVETTE



RASPBERRY
BERET



DEEP
PURPLE



GREEN
MACHINE



BLUE
JEAN



BACK IN
BLACK

**power:
mobility**
The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034

sales@powermobility.com.au

powermobility.com.au

QUANTUM®

#1 FOR REHAB POWER

(TEL) 03 8770 9600 | (FAX) 03 9703 2720

quantumrehab.com

