

**NEW**

# QUANTUM® iLevel® power chairs

New iLevel™ Power Chairs from Quantum Rehab® offer safe eye level operation at the speed of life. iLevel's advanced functionality and stability strengthens independence and quality of life. iLevel improves access to one's environment and enhances social interaction by allowing operation of the power chair with the seat elevated (up to 25.4 cm (10")) while at walking speed (up to 4.83 km/h (3 mph)).



## Introducing a superior seat elevator solution.

- 25.4 cm (10") of elevation while at walking speed (up to 4.83 km/h (3 mph))
- Advanced functionality and stability strengthens daily independence and quality of life for the customer
- Only available on the new Q6 Edge® 2.0 with TRU-Balance® 3 Seating System
- Retrofittable on the Q6 Edge® 2.0 base
- Minimal additional cost



**QUANTUM®**

(Tel) 03 8770 9600 • (Fax) 03 9703 2720

20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

[www.pridemobility.com.au](http://www.pridemobility.com.au)

