


# J 623 2.0



AVAILABLE WITH 2.5 KPH  
12"  
SAFE SEAT ELEVATION



## SOLID PERFORMANCE. SIMPLE CONTROL.

The Jazzy® 623 2.0 is a high-performance mid-wheel drive chair engineered for active users seeking confident, responsive mobility. Powered by robust motors capable of speeds up to 10 km/h (6 mph), it offers smooth handling and strong outdoor performance. Standard TRU-Balance® 3 tilt (50°) provides comfort and positioning, while optional 25.4 cm (12") lift elevation supports greater access and independence.

The VR2 controller enables two powered seating functions, offering simple and reliable controls for everyday driving and power positioning.



**QUANTUM®**  
#1 FOR REHAB POWER



quantumrehab.com.au

## Specifications

<b>DRIVE WHEELS (Mid-Wheel Drive)</b>	35.5 cm (14") Pneumatic or 14x4" wide Flat Free Tyres <sup>2</sup>
<b>CASTER WHEELS</b>	<b>Front:</b> 15.2 cm (6") Solid <b>Rear:</b> 15.2 cm (6") Solid
<b>ANTI -TRIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 10 km/h (2.5 km/h w/ elevation)
<b>GROUND CLEARANCE</b>	6.8 cm (2.7") (battery tray)
<b>TURNING RADIUS</b>	52 cm (20.5")
<b>OVERALL LENGTH</b> <sup>2</sup>	90 cm (35.5") without foot riggings
<b>BASE WIDTH</b>	61.9 cm (24.4")
<b>SEATING SIZES (TRU-Balance® 3)</b>	
<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 12-22" (30.5cm - 55.8cm) <b>Depth:</b> 12-22" (30.5cm - 55.8cm)
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22"(30.5cm - 55.8cm) <b>Depth:</b> 12-22"(30.5cm - 55.8cm)
<b>TRU-Balance® 3 Power Lift &amp; Tilt</b>	<b>Width:</b> 12-22"(30.5cm - 55.8cm) <b>Depth:</b> 12-22"(30.5cm - 55.8cm)
<b>SEAT-TO-FLOOR HEIGHTS ( TRU-BALANCE® 3)</b>	
<b>TRU-Balance® 3 Static Seat</b> <sup>5,7</sup>	16.5-18.5" (42cm - 47cm)
<b>TRU-Balance® 3 Power Tilt</b> <sup>5,7</sup>	17.5 -19.5" (44.5cm - 49.5cm)
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b> <sup>5,7</sup>	17.5-19.5" (44.5cm - 49.5cm)
<b>BATTERY SIZE</b>	MK 50Ah(2)
<b>BATTERY WEIGHT</b> <sup>6</sup>	17.2 kg each
<b>AVAILABLE ELECTRONICS</b>	90A VR2 6-Key
<b>RANGE PER CHARGE</b> <sup>3</sup>	Up to 24.9 km
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-pole motors
<b>WEIGHT CAPACITY</b> <sup>4</sup>	136 kg
<b>BASE WEIGHT</b>	185 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 136 kg weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 1.5kg. based on manufacturer.
- 7) Seat-to-floor heights are measured with the unit unloaded.

Call for availability


## Included Features

- Quick height Adjustable (QHA) armrest
- TB3 Tilt-TRU-Balance 50 degree tilt
- Swingaway legrests with angle adjustable foot plate or center-mount footboard
- Choice of either TRU-Comfort back or upgrade to ADI backs
- Choice of Stealth cushion
- Stealth headrest

## Options

- TRU-Balance® 3 Power Positioning: allows 30.5 cm 12" of seat elevation while at speeds of (2.5 km/h)
- 10cm (14") wide flat-free or pneumatic knobby extended traction tyres<sup>9,4</sup>
- Choice of Stealth cushion
- Stealth headrest
- 8" caster wheels all round

## Colors



**power:mobility**  
The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034  
sales@powermobility.com.au • powermobility.com.au



**QUANTUM®**  
#1 FOR REHAB POWER

 **PROUDLY BUILT IN AUSTRALIA**

(03) 8770 9600 | quantumrehab.com.au