

"FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT'S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE"

Sue
United Kingdom

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.

REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be made-to-order.

UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



Frontier V4 RWD
**ADVENTURE GAINED
NOTHING LOST**



07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034
sales@powermobility.com.au
powermobility.com.au



V6/AUST/NOV21

TAILORED FOR YOUR LIFESTYLE



AUSTRALIAN MADE





FRONTIER V4 REAR WHEEL DRIVE

A powered wheelchair for people who know what they want. The thinking person's chair, the V4 RWD is designed for people wanting to get around the streets or off-road, as well as indoors. Choose off-road or hybrid wheels to ensure you get the best powerchair for your lifestyle.

"WE'VE BEEN
CREATING
MADE-TO-ORDER
WHEELCHAIRS
FOR OVER 25
YEARS, SO THAT
YOU CAN LIVE THE
LIFE YOU
CHOOSE"

Ashley Daff
R&D Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

QUALITY BUILT

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

Basic specifications* (AT wheels and castors):

- Length: 930mm
- Width: 710mm (670mm with Hybrid wheels)
- Seat to floor height: 435mm
- Turning radius: 520mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" black, Hybrid 14"x3" black and grey (solid or pneumatic)
- Castors: 220mm (2.8/2.5-4), black or grey (solid or pneumatic)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:**

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

Base options:



FOREST GREEN - BLACK AT WHEELS & RIMS



SUNRISE YELLOW - BLACK HYBRID WHEELS (PNEUMATIC OR SOLID)



BLAZE BLUE - GREY HYBRID WHEELS (PNEUMATIC OR SOLID)

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Because the world's not flat



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations