

“FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT’S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE”

Sue
 United Kingdom

FULL LIGHTING PACKAGES

Been out having too much fun and it’s getting late? Put your lights on.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



POWER OPTIONS

Consider a seat lift to help you reach, transfer and face life head on, or power tilt for improved manoeuvring over obstacles.

UP CLOSE AND PERSONAL

Side transfers aren’t the only reason you need to flip your armrest out of the way.



Frontier V6 Compact 40
POWER IN A COMPACT PACKAGE

power mobility
 The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034
 sales@powermobility.com.au
 powermobility.com.au



V6/AUST/NOV21

TAILORED FOR YOUR LIFESTYLE



AUSTRALIAN MADE

Frontier
V6™



FRONTIER V6 COMPACT 40

The baby of the family, the Compact 40 is ideal for negotiating tight turns. Small but perfectly formed, the only thing that we've reduced is the size. Still the same frame, motors and seating. We just made it smaller.

"WE'VE BEEN CREATING MADE-TO-ORDER WHEELCHAIRS FOR OVER 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff
R&D Director




ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

QUALITY BUILT

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

Basic specifications:*

- Length: 910mm
- Width: 655mm
- Seat to floor height: 435mm
- Turning radius: 490mm
- Weight limit 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Hybrid 14"x3" grey or black (solid or pneumatic)
- Castors: 200x50 grey and black (solid or pneumatic)
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide)
- Crash test compliant tie down points.

Rehab positioning:

- Power options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Center post, swingaway, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin and environmental controls
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:**

- Every day off-road recreation
- Urban outdoor living and working
- Roomy indoor spaces
- Tight indoor spaces.

Base options:



ULTRA VIOLET - GREY WHEELS (PNEUMATIC OR SOLID)



HOT PINK - GREY WHEELS (PNEUMATIC OR SOLID)



SILVER STREAK - BLACK WHEELS (PNEUMATIC OR SOLID)

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Because the world's not flat



FAMILY OF BRANDS

Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations