

### WHEELCHAIR RANGE

aspire for...

Comfort
Posture
Mobility













EVOKE 2

REHAB RX

**ASSIST 2** 

NSIT 2 LITE TRANSIT



www.aspirecare.com.au 1300 133 120



### aspire EVOKE 2 Lightweight & Highly Adjustable

The Aspire EVOKE 2 and EVOKE 2 HD are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for every user. A comprehensive range of compatible accessories is available for increased independence and enhanced clinical outcomes.



#### **FEATURES**

Frame: Folding Lightweight Aluminium

Legrest: Swing-away & Removable

Height & Angle Adjustable **Footrest:** 

Armrest: Swing-back & Removable

Quick Release, Solid 24" Tyres, Wheels:

Height & Depth Adjustable

**Castors:** 8", Height & Angle Adjustable

**Brakes:** Easy Access Push-to-Lock Brakes

Depth & Height Adjustable Seat:

**Backrest:** Height, Tension, Angle & Depth Adjustable

Dacron Fire Retardant Fabric **Upholstery:** 

Contoured Overlay is optional, **Cushion:** 

Standard cushion Included

**OPTIONAL ACCESSORIES** 

Elevating Legrest

**Elevating Legrest** 

Stump Supports

Advanced Articulating

Oxygen Bottle Holder

Weight: From 14.8kg (SWL 140kg - HD 180kg)




**HEIGHT ADJUSTABLE** 

**ADJUSTABLE** 

#### OPTIONAL CONTOURED PADDED OVERLAY

This offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. The contouring assists in positioning and pressure distribution, while the fabric allows air movement, increasing breathability. It is easily removed for washing.



#### **CUSTOMISABLE BACKREST**

Tension Adjustable

✓ Angle Adjustable

Height Adjustable
















SEAT WIDTH	CODE
350 mm (14")	MWS449500
400 mm (16")	MWS449510
450 mm (18")	MWS449520
500 mm (20")	MWS449530
570 mm (22")	MWS449540 (HEAVY DUTY model)
610 mm (24")	MWS449550 (HEAVY DUTY model)

#### \*For Detailed Features, Specifications and Accessories see pages 14 - 20 Note: Side transfers are not possible when gutter armrests are attached

Brake Extensions

Drum Brakes

Lateral Supports

Gutter Arm Rests\*

Contoured Padded Overlay

Aspire EVOKE 2 is crash tested to AS/NZS3696.19

Anti-tip Bars

Safety Belt

Camber Kit

Calf Straps

IV Pole

Headrest Assembly

18" - 21"

Extended Back Canes

**ADJUSTABLE** 

# aspire EVOKE 2 JNR




2 have been reduced in size to suit the unique needs of children, including smaller 16 inch rear wheels.

#### **CUSTOMISABLE BACKREST**

- Depth Adjustable
- (V) Angle Adjustable
- (V) Height Adjustable









HEIGHT ADJUSTABLE DEPTH ADJUSTABLE

Frame:

Folding Lightweight Aluminium

Legrest:

Swing-away & Removable

Footrest:

Height & Angle Adjustable

Armrest:

Height Adjustable &

Wheels:

Quick Release, Solid 20" Tyres

**Castors:** 

8", Height & Angle

Adjustable

Easy Access Push-to-lock

**Brakes:** Seat:

Height & Depth Adjustable

**Backrest:** 

Height, Angle & Depth

Adjustable

**Upholstery:** 

Contoured Overlay is Optional, Standard Cushion Included

Cushion: Included Weight: From 14kg

\*For Detailed Features, Specifications and Accessories see pages 14 - 20











	SEAT WIDTH	
$\leq$	300 mm (12")	
	A : FVOVE 2 IND :	

Aspire EVOKE 2 JNR is crash tested to AS/NZS3696.19

## aspire REHAB RX JNR

The Aspire REHAB RX JNR is designed for advanced postural support and pressure care for children via easy repositioning and advanced adjustability.

ADVANCED COMFORT Breathable PU coated multi-stretch fabric on seat and back cushions ensure optimal skin microclimate, comfort and improve pressure care outcomes.

TILT-IN-SPACE (TIS) Functionality promotes improved positioning. As well as this, it allows body

areas.

CUSHIONING Multi-layer, castellated seat and back cushions with internal postural reinforcement optimise comfort and positioning.







**FEATURES** 

Frame:

Legrest:


















TILT-IN-SPACE







Height & Angle Adjustable Footrest:

Rigid Steel

Removable, Height & Depth Adjustable **Armrest:** 

Quick Release, Solid 20" Tyres, Wheels: Height & Depth Adjustable

**Castors:** 6", Height Adjustable

Attendant Operated Drum Brakes **Brakes:** & Push-to-lock Handbrakes

Depth & Width Adjustable Seat:

**Backrest:** Height, Tension & Angle Adjustable

Height, Depth, Angle & Width **Headrest:** Adjustable, Removable

PU Breathable Multi-stretch **Upholstery:** 

Fabric Covers

Multilayer Castellated Back & **Cushion:** Seat Cushions in Breathable PU

Weight: From 39.5kg



CODE

Aspire REHAB RX JNR is crash tested to AS/NZS3696.19















## aspire REHAB RX Advanced Tilt-In-Space

The REHAB RX is geared towards the high needs, rehabilitation environment. Highly configurable and with pressure care cushioning surfaces, the REHAB RX is designed for the user with more complex



#### **FEATURES**

**Backrest:** 

Frame: Rigid Steel Frame

Swing-away, Removable, Padded Legrest:

& Elevating

**Footrest:** Height & Angle Adjustable

Removable, Height & Depth Adjustable **Armrest:** 

Quick Release, Solid 24" Tyres. Wheels: Height & Depth Adjustable

6", Height Adjustable **Castors:** 

Attendant Operated Drum Brakes **Brakes:** 

& Push-to-lock Handbrakes

Height, Tension & Angle Adjustable

Depth. Width & Height Adjustable Seat:

Height, Depth, Angle & Width

**Headrest:** Adjustable, Removable

PU Breathable Multi-Stretch Fabric Covers **Upholstery:** 

Multilaver Castallated Back and Seat Cushion **Cushion:** 

Covered in Breathable PU Fabric

Weight: From 41kg

\*For Detailed Features, Specifications and Accessories see pages 14 - 20

#### **OPTIONAL ACCESSORIES**

A wide range of accessories to meet clinical needs

- Base Frame Extension
- 20" or 16" Alternative
- Wheels
- Therapist Adjustment Kit
- Stump Supports
- Safety Belt
- · Oxygen Bottle Holder
- Meal Tray
- Anti-tip Bars
- IV pole

### **HEAD SUPPORT** Easy-to-adjust neadrest offers the user superior support & comfort. Breathable hygiene PU surfaces, easy to clean Swing-back armrests • CUSHION (height and depth adjustable option available) Quick release wheels, height & depth removable and padded Height & angle Height adjustable 6" adjustable footrests

#### **SUPPORT SURFACES:**

- Multi-layer foam construction with castellated top layer for enhanced comfort, pressure redistribution and durability
- Breathable PU multi-stretch cover for improved pressure redistribution and skin micro-climate management

#### **BACKREST - CONTOURED:**

- Gentle lumbar build up for postural support
- Moderate lateral supports for assisted midline sitting

#### **CUSHION:**

- Integrated positioning skeleton for enhanced sitting posture
- Fabric tensioners provide a smooth support surface
- Gentle abductor and adductor contouring provide improved postural support

#### ADJUSTABLE ELEVATING LEGRESTS

Elevating legrests with height & angle adjustable footrests.

#### INFECTION CONTROL Easy to clean surfaces.

**CUSTOM COMFORT** After-market backrests

& cushions ensure flexibility.



#### **SAFETY IN DESIGN**

Padded fittings and supports keep the user safe and secure at all times.

Aspire REHAB RX is crash tested to AS/NZS3696.19







SEAT WIDTH	CODE
410 - 460 mm (16" - 18")	MWS449710
460 - 510 mm (18" - 20")	MWS449720
510 - 560 mm (20" - 22")	MWS449740










#### THE HIGHLY ADJUSTABLE FEATURES FULFILL MANY CLINICAL REQUIREMENTS GEARED TOWARDS COMFORT AND SUPPORT





castors



adjustable

## aspire REHAB RS Classic Tilt-In-Space

The REHAB RS is an ideal solution for Aged Care and Community users who requires increased postural support and bodyweight redistribution. A reclined position is easily achieved, making the REHAB RS also a great option for those who fatigue easily.



#### **FEATURES**

Rigid Steel Frame Frame:

Swing-away, Removable, Padded & Elevating Legrest:

**Footrest:** Height & Angle Adjustable

Removable, Height & Depth Adjustable Armrest:

Quick Release, Solid 24" Tyres, Wheels:

Height & Depth Adjustable

**Castors:** 6", Height Adjustable

Attendant Operated Drum Brakes **Brakes:** 

& Push-to-lock Handbrakes

Seat: Depth, Width & Height Adjustable

Height, Tension & Angle Adjustable **Backrest:** 

Height, Depth, Angle & Width **Headrest:** Adjustable, Removable

**Upholstery:** PU Breathable Multi-Stretch Fabric Covers

**Cushion:** Standard Foam Comfort Cushions on Back & Seat

Weight: From 40.7kg

\*For Detailed Features, Specifications and Accessories see pages 14 - 20

#### **OPTIONAL ACCESSORIES**

A wide range of accessories to meet clinical needs

- Base Frame Extension
- 20" or 16" Alternative
- Wheels
- Therapist Adjustment Kit
- Stump Supports
- Safety Belt
- Oxygen Bottle Holder
- Meal Tray
- Anti-tip Bars
- IV pole



#### **INCREASED HEAD SUPPORT**

The easy-to-adjust, supportive headrest offers the user superior head support, while the gently contoured backrest facilitates posture and positioning.

#### **SUPPORTIVE & COMFORTABLE**

Padded armrests, elevating legrests & a supportive headrest ensure that users of the REHAB RS are supported and comfortable.

#### **ADVANCED COMFORT**

Breathable PU coated multi-stretch fabric on seat and back cushions ensure optimal skin microclimate, comfort and improve pressure care outcomes.



#### **EASY TO CLEAN**

Wipe clean surfaces ensure infection control compliance.



#### **SAFETY AT ITS CORE**

Padded armrests. elevating legrests & a supportive headrest ensure that users of the REHAB RS remain safe and secure.



Aspire REHAB RS is crash tested to AS/NZS3696.19







SEAT WIDTH
410 - 460 mm (16" - 18")
460 - 510 mm (18" - 20")

510 - 560 mm (20" - 22")

MWS449600 MWS449610 MWS449620

CODE





DVA Contracted



#### THE HIGHLY ADJUSTABLE FEATURES FULFILL MANY CLINICAL REQUIREMENTS GEARED TOWARDS COMFORT AND SUPPORT











#### LOW TO MEDIUM LEVEL REHABILITATION SOLUTION

Tilt-In-Space (TIS) functionality promotes increased support and improved positioning. This allows caregivers to position the user for optimal comfort. Adjustable, elevating legrests allow for precise leg positioning.



1300 133 120

### aspire ASSIST 2 Everyday

The Aspire ASSIST 2 is a versatile and robust self-propelling wheelchair, designed to meet the rigorous demands of Community, Hire, Hospital and Aged Care environments. The Aspire ASSIST 2 is a reliable and comfortable solution with various accessories available to suit the needs of individuals and caregivers alike.





#### **FEATURES**

Footrest:

**Castors:** 

Folding Steel Frame Frame: Swing-away, Removable, Legrest:

(Elevating optional accessory) Height Adjustable

**Armrest:** Swing-back, Removable

Wheels: Quick Release, Solid 24" Tyres

8", Height Adjustable

**Brakes:** Push-to-lock

Fixed, Light Padding Seat: **Backrest:** Fixed, Light Padding

**Upholstery**: Dacron Fire Retardant Fabric

**Cushion:** Optional Accessory

Weight: From 17.3kg

#### **ROBUST & RELIABLE**

With a strong steel frame, removable swing-away legrests and quick release self-propelling rear wheels, the Aspire ASSIST 2 is a practical yet robust patient transport solution.

#### **CATERS TO A WIDE VARIETY OF ENVIRONMENTS**

The large, puncture-proof wheels ensure that the Aspire ASSIST 2 handles varied terrain with ease.

#### LOW MAINTENANCE FOLDING

Solid rear tyres and a centre cross brace design offer a low maintenance and portable mobility solution within a facility, hire or community environment.

SEAT WIDTH	CODE
350 mm (14")	MWS449300
400 mm (16")	MWS449310
450 mm (18")	MWS449320
500 mm (20")	MWS449330









\*For Detailed Features, Specifications and Accessories see pages 14 - 20



The Aspire ASSIST 2 WARD Vinyl wheelchair has all the features of the ASSIST 2, but with durable hospital grade vinyl, suitable for the hospital & clinic environment. The ASSIST 2 WARD Vinyl has been carefully designed with fully welded seams and concealed mounting of upholstery to ensure that rigorous infection control standards can be met.





#### **FEATURES**

Frame: Folding Steel Frame Swing-away, Removable Legrest: (Elevating optional accessory)

Height Adjustable Footrest:

Armrest: Swing-back, Removable

Wheels: Quick Release, Solid 24" Tyres

8", Height Adjustable **Castors:** 

Push-to-lock **Brakes:** 

Fixed, Padded, Fire Retardant Vinyl Seat:

**Upholstery:** Hospital Grade Fire Retardant Vinyl

Fixed, Padded, Fire Retardant Vinyl

**Cushion:** Optional Accessory

Weight: From 17.3kg

Simple Frame Fix Mod to Stop

Chair Being Folded (Available as option on request)

#### **INFECTION CONTROL FOCUS**

Durable hospital vinyl that is easy to wipe clean

Fully welded seams to inhibit bacterial growth


Concealed mounting of upholstery to prevent fluid ingress

400	/	~ - /	~ I II	<b>-</b> 11	- > /
ASS	⊢ I '	$\sim$		01	ıv
~~		ューヽ	_		

The centre cross-brace of the ASSIST 2 WARD can be customised to fix the frame for Hospital Asset Security.

#### **DURABLE TRANSPORT SOLUTION**

The ASSIST 2 WARD Vinyl is a versatile transport wheelchair. Quick release wheels with solid, maintenance-free tyres provide on-call reliability and are suitable for use in common areas, corridors, wards, bathrooms and outdoors.







1300 133 120











# aspireTRANSIT 2 Transporter

The Aspire TRANSIT 2 wheelchair is a versatile and robust transport wheelchair that is propelled by an attendant or carer. The Aspire TRANSIT 2 is a reliable and comfortable solution, with various accessories available to suit the needs of individuals and carers in the Community, Hospital or Aged Care environments.





#### **FEATURES**

**Brakes:** 

Frame: Folding Steel Frame Swing-away, Removable, Legrest: (Elevating optional accessory)

Height Adjustable Footrest:

**Armrest:** Swing-away, Removable

Wheels: Quick Release, Solid 12" Tyres

Attendant Push-to-lock

**Castors:** 8", Height Adjustable

Seat: Fixed, Light Padding

**Backrest:** Fixed, Light Padding

Dacron Fire Retardant Fabric **Upholstery:** 

**Cushion:** Optional Accessory

Weight: From 16.4kg

#### **SWING-AWAY ARMRESTS**

The padded armrests are swing-away and removable to facilitate the user when transferring in and out of the chair.

#### **DURABLE TRANSPORT SOLUTION**

The Aspire TRANSIT 2 is a versatile transport wheelchair, with rear 12 inch puncture-proof, transit wheels to allow an attendant or carer to easily navigate a client to their destination.

#### **FOLDING DESIGN**

The centre cross-brace design allows for the Aspire TRANSIT 2 to be folded up for unoccupied transport or storage within a facility, hire or community environment.

SEAT WIDTH	CODE
400 mm (16")	MWS449410
450 mm (18")	MWS449420
500 mm (20")	MWS449430











\*For Detailed Features, Specifications and Accessories see pages 14 - 20

### aspire LITE TRANSIT

The Aspire LITE TRANSIT is an lightweight wheelchair ideal for carer assisted transport. The compact lightweight frame makes it easy to transfer in and out of the car. Attendant brakes and 12 inch rear wheels assist the carer in negotiating indoor and outdoor terrain smoothly and safely. While the height adjustable armrests allow the user close access to tables or quick adjustment for comfort.



The Aspire LITE is an 'ultra lightweight' wheelchair that is ideal for occasional use. The Aspire LITE is designed with the carer in mind. The compact, lightweight frame makes it easy to transfer in and out of a car. The push button, height adjustable armrests make it versatile, allowing the user to be positioned at a table and quickly adjusted for comfort.



#### **FEATURES**

**SAFETY** 

Folding Aluminium Frame Frame: Legrest: Swing-away & Removable

**Armrest:** Height Adjustable Wheels: Solid 12" Tyres **Castors:** 6" Height

Attendant Brakes **Brakes:** & Push-to-lock Brakes

**Upholstery:** Padded Dacron Fabric Weight: Chair Weight 11.6kg

### S.W.L

HEIGHT

ADJUSTABI F

ARMRESTS &

#### **FEATURES**

Frame: Folding Aluminium Frame Legrest: Swing-away & Removable

Height Adjustable **Armrest:** Wheels: Solid 8" Tyres **Castors:** 6" Height

Push-to-lock **Brakes:** 

**SAFETY & STORAGE** 

Safety lap belt & rear

**ULTRA LIGHTWEIGHT** 

minimal manual handling.

At 9.5kg, allows for

storage pocket

Upholstery: Padded Dacron Fabric Chair Weight 9.5kg Weight:



HEIGHT

ADJUSTABLE

#### **TRANSPORTABLE**

Safety lap belt & Easily folded due to it's attendant brakes centre cross brace design. The attendant handles can LIGHTWEIGHT also be folded down, taking

Reduced manual up minimal car space. handling effort.

CHAIR	SEAT WIDTH	CODE
TRANSIT LITE	460 mm (18")	MWS442301
LITE	460 mm (18")	MWS442300











up minimal car space.

**TRANSPORTABLE** 

Easily folded due to it's

centre cross brace design.

The attendant handles can

also be folded down, taking





					SUITS				
DESCRIPTION									B.
	Evoke 2	Evoke 2 JNR	Rehab RX	Rehab RX	Rehab RS	Assist 2	Transit 2	Lite Transit	Lite
Aluminium	•	JNR	JNR						
Steel									
Folding									
SEAT									
Depth adjustable									
Width adjustable									
Cushion included									
WHEELS									
Quick release	•	•							
Solid tyres									
24" wheels									
20" wheels									
12" wheels									
8" wheels									
LEGRESTS									
Swing-away									
Removable									
Padded									
Elevating									
ARMRESTS									
Swing-away									
Removable									
Height adjustable									
Depth adjustable									
BACKREST									
Height adjustable back canes									
Tension adjustable back caries									
BRAKES									
Drum Brakes									
Attendant Brakes									
OTHER FEATURES Self-propelled									
Attendant-propelled only									
Tilts									
Reclines									
Crash tested									

### TECHNICAL SPECIFICATIONS

CODE	SEAT WIDTH	OVERALL WIDTH	SEAT DEPTH	SEAT HEIGHT	BACKREST HEIGHT	ARMREST HEIGHT	MINIMUM CHAIR WEIGHT	S.W.L
EVOKE 2				<b>\$</b>	<b>†</b>	<b>†</b>	KG	KG
MWS449500	350 mm (14")	560 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 14.8 kg	140 kg
MWS449510	400 mm (16")	610 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 14.9 kg	140 kg
MWS449520	450 mm (18")	660 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 15.0 kg	140 kg
MWS449530	500 mm (20")	710 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 15.1 kg	140 kg
MWS449540 (HD)	570 mm (22")	780 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 16.0 kg	180 kg
MWS449550 (HD)	610 mm (24")	820 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 16.2 kg	180 kg
EVOKE 2 JNR				<b>\$</b>	<b>†</b>	<b>†</b>	KG	KG
MWS449560	300 mm (12")	510 mm	300 - 375 mm	360 - 410 mm	430 - 490 mm	205 - 225 mm	14 kg	140 kg
REHAB RX JNR	<b>≟</b> <b>⊟</b>		<b>*</b>	<b>\$</b>	<b>*</b>	<b>†</b>	KG	KG
MWS449630	320 - 370 mm (13"- 15")	600 mm	440 - 490 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	39.5 kg	180 kg
REHAB RX				<b>↑</b>	<b>\$</b>	<b>\$</b>	KG	KG
MWS449710	410 - 460 mm (16" - 18")	690 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	41 kg	180 kg
MWS449720	460 - 510 mm (18" - 20")	740 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	43 kg	180 kg
MWS449740	510 - 560 mm (20" - 22")	790 mm	470 - 560 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	47.2 kg	180 kg
REHAB RS			**	<b>\$</b>	<b>†</b>	<b>†</b>	KG	KG
MWS449600	410 - 460 mm (16" - 18")	690 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	40.7 kg	180 kg
MWS449610	460 - 510 mm (18" - 20")	740 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	42.5 kg	180 kg
MWS449620	510 - 560 mm (20" - 22")	790 mm	470 - 560 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	46.7 kg	180 kg

e: Seat height for REHAB RS, RX, RX JNR, EVOKE 2 & EVOKE 2 JNR is without seat cushion. Standard seat cushion height is 100 mi

www.aspirecare.com.au 1300 133 120

### TECHNICAL SPECIFICATIONS

CODE	SEAT WIDTH	OVERALL WIDTH	SEAT DEPTH	SEAT HEIGHT	BACKREST HEIGHT	ARMREST HEIGHT	MINIMUM CHAIR WEIGHT	S.W.L
ASSIST 2	<b>→</b>	<b>→</b>		<b>\$</b>	<b>†</b>	<b>†</b>	KG	KG
MWS449300	350 mm (14")	570 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.3 kg	140 kg
MWS449310	400 mm (16")	620 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.5 kg	140 kg
MWS449320	450 mm (18")	670 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.7 kg	140 kg
MWS449330	500 mm (20")	720 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.9 kg	140 kg
MWS449340 (DELUXE VINYL)	350 mm (14")	570 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.4 kg	140 kg
MWS449350 (DELUXE VINYL)	400 mm (16")	620 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.6 kg	140 kg
MWS449360 (DELUXE VINYL)	450 mm (18")	670 mm	400 mm	460 - 500 mm	440 mm 210 mm		17.8 kg	140 kg
MWS449370 (DELUXE VINYL)	500 mm (20")	720 mm	400 mm	460 - 500 mm	440 mm 210 mm		18 kg	140 kg
MWS449380 (HD DELUXE VINYL)	570 mm (22")	790 mm	400 mm	460 - 500 mm	440 mm	210 mm	18.9 kg	180 kg
MWS449390 (HD DELUXE VINYL)	610 mm (24")	830 mm	400 mm	460 - 500 mm	440 mm	210 mm	19.2 kg	180 kg
TRANSIT 2	<b>*</b>	<b>*</b>	**	<b>†</b>	<b>†</b>	<b>†</b>	KG	KG
MWS449410	400 mm (16")	540 mm	400 mm	460 - 500 mm	440 mm	210 mm	16.4 kg	140 kg
MWS449420	450 mm (18")	590 mm	400 mm	460 - 500 mm	440 mm	210 mm	16.6 kg	140 kg
MWS449430	500 mm (20")	640 mm	400 mm	460 - 500 mm	440 mm	210 mm	16.8 kg	140 kg
LITE TRANSIT	<b>→</b>			<b>\$</b>	<b>†</b>	<b>†</b>	KG	KG
MWS442301	460 mm (18")	600 mm	400 mm	460 mm	465 mm	235 mm	11.6 kg	120 kg
LITE				\$	<b>‡</b>	<b>†</b>	KG	KG
MWS442300	460 mm (18")	550 mm	400 mm	460 mm	465 mm	235 mm	9.5 kg	120 kg

### **ACCESSORIES**









BASE FRAME EXTENSION KITS

STUMP SUPPORT

**ELEVATING LEGRESTS** 

















**OVERLAY** 

ARMRESTS

CALF STRAP

**GUTTER ARMRESTS** 









SAFETY LAP BELT

**BRAKE EXTENSION LEVERS** 

LATERAL CHEST SUPPORT

TOOL-FREE ADJUSTMENT KIT



**IV POLE** 





**EXTENDED BACK CANES** 

### ACCESSORIES MATRIX

						SUITS			
CODE	DESCRIPTION		Evoke 2	Evoke 2 JNR	Rehab RX	Rehab RX	Rehab RS	Assist 2	Transit 2
ELEVATING L	EGRESTS								
WCA848100	Elevating Legrest - Left - 300/320/350/400/410mm	A			0	0	0	•	
WCA848110	Elevating Legrest - Right - 300/320/350/400/410mm			•	0	0	0	•	
WCA848120	Elevating Legrest - Left - 450/460/500/510mm					0	0		
WCA848130	Elevating Legrest - Right - 450/460/500/510mm		•			0	0	•	
WCA848140	Elevating Legrest - Left - 570/610mm								
WCA848150	Elevating Legrest - Right - 570/610mm		•					•	
STUMP SUPPO	DRTS								
WCA848220	Stump Support - Left - Silver - 300/320/350/400/450mm		•	•	•	•	•	•	•
WCA848230	Stump Support - Right - Silver - 300/320/350/400/450mm		•	•	•	•	•		
WCA848240	Stump Support - Left - Silver - 500/550/570/610mm		•			•	•	•	•
WCA848250	Stump Support - Right - Silver - 500/550/570/610mm	14	•			•	•	•	
WCA848260	Stump Support - Extension Bracket - Silver		•	•				•	•
OXYGEN BOT	TLE HOLDERS								
WCA848290	Oxygen Bottle Holder	1	•	•				•	
WAA695600	Oxygen Bottle Bag				•	•	•		
BRAKE EXTER	NSIONS								
WCA848310	Brake Extension Lever - Left		•	•				•	
WCA848320	Brake Extension Lever - Right	7	•	•				•	•
LATERAL SUF	PPORTS								
WCA848430	Lateral Chest Support - Left	4	•	•					
WCA848440	Lateral Chest Support - Right	No.							
WCA848450	Lateral Chest Support - Left								
WCA848460	Lateral Chest Support - Right								

### ACCESSORIES MATRIX

		SUITS						
CODE	DESCRIPTION	Evoke 2	Evoke 2 JNR	Rehab RX JNR	Rehab RX	Rehab RS	Assist 2	Transit 2
CONTOURED PADDED OVERLAYS								
WCA848330	Contoured Padded Overlay - 350mm Wide							
WCA848340	Contoured Padded Overlay - 400mm Wide	•						
WCA848350	Contoured Padded Overlay - 450mm Wide	•						
WCA848360	Contoured Padded Overlay - 500mm Wide							
WCA848370	Contoured Padded Overlay - 570mm Wide	•						
WCA848380	Contoured Padded Overlay - 610mm Wide						•	
DEPTH & HEIGHT ADJUSTABLE ARMRESTS								
WCA848390	Depth & Height Adjustable Armrest - Left	•						
WCA848400	Depth & Height Adjustable Armrest - Right	•						
GUTTER ARMRESTS *								
WCA848410	Gutter Armrest Pad - Left	•						
WCA848420	Gutter Armrest Pad - Right	•						
DRUM BRAKES								
WCA848510	Rear Wheels w/Drum Brakes Kit - 16in		•	•	•	•		
WCA848520	Rear Wheels w/Drum Brakes Kit - 20in		•	0	•	•		
WCA848500	Rear Wheels w/Drum Brakes Kit - 24in				0	0		
DEPTH EXTENSIONS								
WCA848530	Base Frame Extension Kit - 50mm - Black				•	•		
WCA848540	Base Frame Extension Kit - 100mm - Black				•	•		
OTHER								
WCA848550	Therapist-Friendly Tool-Free Adjustment Kit			•	•	•		
WCA848560	Meal Tray - 320mm to 410mm			•	•	•		
WCA848570	Meal Tray - 460mm to 510mm					•		
WCA848270	Anti-Tip Bars	•	•	0	0	0	•	
WCA848280	IV Pole		•	•		•		
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle	•	•	•	•	•		
WCA848470	Extended Back Cane - Pair - 19/21in		•					
WCA848480	Headrest Assembly - Deluxe Style	•	•					
WCA848580	Calf Strap	•	•	•	•	•		
WCA848490	Camber Kit	•	•					
AAL848590	Security Mod - ASSIST 2 WARD Only							

O - Standard Feature Included \*Note: Side transfers are not possible when gutter armrests are attached




### aspire Associated Products



































07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034 sales@powermobility.com.au • powermobility.com.au

