

Avantgarde XXL 2

The Avantgarde XXL 2 was developed to meet the needs of tall, heavier people, and its strengths include outstanding stability combined with easy operation. The robust base structure offers high reliability in any situation. The choice of high-quality materials keeps the net weight as low as possible.

- User weight up to 180 kg
- Reinforced aluminium tube frame
- Reinforced hollow rims
- Segmented, angle-adjustable foot support with aluminium foot plates



07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034 sales@powermobility.com.au powermobility.com.au

Technical data

Seat width	460-620 mm
Front seat height	440-530 mm
Rear seat height	410-530 mm
Seat depth	400-560 mm
Overall length	790-1,140 mm
Overall width	695-940 mm
Back support height	250-500 mm
Lower leg length	320-550 mm
Approx. weight	17 kg
Transport weights Frame approx. Drive wheel from Foot supports approx.	11.3 kg 2.15 kg 1.2 kg
Max. load	180 kg