

Duro Stroller

Whether navigating city footpaths or the outdoors terrain, Duro is a stable stroller solution with a great range of support and positioning accessories that will enable you to embrace the outdoors.

Configurable rehab stroller for everyday challenges

Extremely agile with **great manoeuvring** due to compact mobility base

Adjustable back angle

Select required **seat angle**

Adjustable knee angle



Generous **growth capacity** provides long-term use for growing users

Individually customisable to suit requirements, with proven accessories from the **Bingo Evo Stroller**

Very stable stroller option



Super compact folding makes the Duro easy to transport and stow in the boot of your car!

Reclining back support means you can relax back in the Duro when being out and about gets tiring!

Colour Options: *(Upholstery must be selected in addition to Stroller)*

Turquoise

Size 1:
2780-8801-339
Size 2:
2780-8802-339




Red
Size 1:
2780-8801-345
Size 2:
2780-8802-345

ACCESSORIES

Standard package includes: hip supports, swivel front wheels, flip-up foot support, soft accessory bag, tie-down kit and retainer for accessories.

Head Support Pads
1: 2880-4461-000
2: 2880-4462-000

Head Support Speakers
2880-8330-000
(Head support pads required)



Occipital Head Support Spacer: 2010-4101-099
Wipe Down: 2010-4101-209

Std Lateral Head Support Spacer: 2010-4201-099
Wipe Down: 2010-4201-209




Std Lateral Support
1: 2880-3010-000
2: 2880-3020-000
3: 2880-3030-000

Contoured Lateral Support: 2006-3030-099
Teardrop Lateral Support: 2006-3110-099



Spex Centrepoint Harness
XS: 3780-4301-209
S: 3780-4301-009

5-Point Harness:
3780-4101-209



Ankle Foot Fast
1: 2880-6101-000
2: 2880-6102-000

Therapy Tray
1: 2880-7111-000
2-3: 2880-7112-000
Infill: 2880-7180-000



Full range of accessories available on page 17

Size	Code	Weight Limit	Seat Width	Seat Depth	Seat Angle	Backrest Height	Backrest Angle	Lower Leg Length	Wheel Diameter	Folded Size (L x W x H)
1	2780-0001-052	75 kg	23 - 34 cm	26 - 40 cm	16° or 20°	45 - 63 cm	78 - 120°	20 - 44 cm	20 cm front, 30 cm rear	84x63x32 cm
2	2780-0002-052	75 kg	28 - 39 cm	31 - 45 cm	16° or 20°	62 - 83 cm	78 - 120°	20 - 44 cm	20 cm front, 30 cm rear	84x68x38 cm