

SOLAX CHARGE

AUTOMATIC FOLDING MOBILITY SCOOTER



**ADJUSTABLE
REAR SUSPENSION**



**REMOTE CONTROL
AUTOMATIC FOLDING**



**USB PORT FOR
CHARGING PHONE**



SOLAX CHARGE

AUTOMATIC FOLDING MOBILITY SCOOTER

The Solax Charge is the first Solax that features a USB charging port that allows you to charge your phone or iPad. Other tech features include automatic folding, LED headlights, rear suspension, a digital dashboard that displays your speed, trip distance and battery life. A re-designed seat boasts a higher backrest with extra padding for additional comfort, while the armrests fold up for easier access.

SPECIFICATIONS

Weight With Batteries	26kg	Tyres Solid - Puncture Proof
Battery Weight	1.7kg	Dimensions LxWxH (mm) 1000x550x865
Heaviest Part	24.3kg	Folded Dimensions (mm) 410x550x690
Weight Capacity	125kg	Distance Between Armrests 480mm
Adjustable Tiller Angle	Yes	Drive Wheel 170mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance 38mm
Battery	24V/10ah Lithium	Seat Dimensions 430mm x 340mm
Battery Range	Up To 15km	Seat Height From Ground 530mm
Maximum Speed	6km/h	Slope Ability 12°
Motor	120W	Turning Circle 1250mm
Number of Components	1 + Battery	Suspension Rear Adjustable



07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034
sales@powermobility.com.au • powermobility.com.au

Colours Subject to Availability

