

SOLAX HOIST

ELECTRIC POWERED HOIST



LIFTS UP TO 30KG



SAFELY LIFT YOUR SCOOTER INTO YOUR CAR



LIGHTWEIGHT COMPACT DESIGN



SOLAX HOIST

ELECTRIC POWERED HOIST

The Solax Portable Electric Hoist helps safely lift your travel scooter in and out of the boot of your car. Made from light telescopic aluminium, this hoist weighs only 13kg which means it can be easily lifted in and out of your vehicle. The hoist is powered by your scooter's battery, making this accessory the perfect companion for your Solax. An external battery pack is also available.

SPECIFICATIONS

Weight Capacity30kg

Total Weight13kg

Dimensions LxWxH (mm)700x850x1780

Maximum SlopeUse on level ground

Lifting Height1500mm

Dimensions When Folded (mm)230x400x700

BatteryPlugs into scooter's battery

Safety MechanismLockable castors

power mobility
The power to move you

07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034

sales@powermobility.com.au • powermobility.com.au