



LIFECOMFORT[®] FOAM MATTRESSES

Foam Support Surfaces for Community, Aged Care and Acute Care



ABOUT LIFECOMFORT® MATTRESSES

LIFECOMFORT® MATTRESSES have set the benchmark in sleep surface quality and performance in Australia for more than a decade.

ALL LIFECOMFORT® MATTRESSES are fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care, Aged Care and up to Acute Care environments.




Multiple foam layering combinations provide maximum comfort, support and pressure redistribution across a wide range of user shapes and sizes.

HOSPITAL
GRADE COVERS
ON ALL
LIFECOMFORT®
MATTRESSES



STANDARD INCLUSIONS

- **Durable Hospital Grade Foam** that is treated for antimicrobial and antibacterial properties.
- **Breathable, 4-Way Stretch Waterproof Cover** to reduce shear, friction and maintain skin microclimate.
- **Medical Grade Waterfall Covers** are fully-welded over the zipper to maximise infection control.
- **Asset Management Ready** with unique batch and serial numbers printed on all covers



LIFE CYCLE & COMPLIANCE


- **Warranty:** 5 yrs Foam | 4 yrs Cover
- **Machine Washable Covers:** Up to 95°
- **Safe Working Load Verified:** AS/NZ 2282
- **Durability + Performance:** ISO 3385
- **Flame Retardant Foam + Covers:** AS/NZ 4088.1 Cribs 5.
- **TGA:** Lifecomfort Mattresses are registered in the ARTG, AUST R ARTG 288486.

CUSTOMISATION

- **Unique User Needs:**
Palliative, unusual body shapes, zoning needs, Roho inserts or specific mattress sizes.
- **Facility Asset Management:**
Customised labelling for institutional orders.

PRESSURE CARE RISK LEVELS*

LIFECOMFORT® MATTRESSES cater to a broad spectrum of Pressure Care risk ratings from Low to High, as illustrated below:




*Pressure care rating references provide general guidelines only, patient specific assessment and correct nursing care are essential for maximising pressure care outcomes.

COMMUNITY AND AGED CARE


COMFORT AND PRESSURE REDISTRIBUTION

Preventative mattress solutions for low to medium pressure relief and maximum comfort. Foams and covers are durability tested to meet the needs of community and modern aged care environments. All mattresses in this range are engineered to move with adjustable beds. Safe Patient Weight ranges have been designed to suit typical Community and Aged Care users maximising sleep comfort and immersion.




1. IMMERSIVE SURFACES

- Envelops vulnerable and bony prominences.
- Evenly distributes surface contact for shear mitigation and improved skin microclimate.
- Contours and cradles around the user for personalised pressure relief and comfort.




2. COMFORT LAYERS

- Provides an additional layer of immersive comfort and envelopment.
- Maintains postural integrity when profiling with an adjustable bed.



3. SUPPORT LAYERS


- Supportive barrier for enhanced comfort in the seated/fowler position.
- Facilitates top layer compression for improved pressure relief.



4. INTEGRATED SIDE WALLS*

- Improves bed mobility offering a firm base of support for entry and exit.
- Side wall firmness increases bed edge awareness for improved falls prevention.

CLINICAL BENEFITS



*LIFECOMFORT Pressure Reducing and Premier Care


LIFECOMFORT® STANDARD

This **LIFECOMFORT® STANDARD** dual layered mattress with visco-elastic memory foam upper, is an economical comfort and care solution for in home use and lighter weight users. The soft and supportive patient surface is ideal for the implementation of palliative care in the home.

COMMUNITY

LOW - MEDIUM

AVAILABLE WITH SAFETY SIDES
An adaptation of the LIFECOMFORT Standard Mattress to include side wall support to help minimise the risk of falls from bed.



FOAM LAYERS

- IMMERSIVE SURFACE
Immersive visco-elastic
- COMFORT, SAFETY & SUPPORT BASE
Supportive high resilience foam


BEM046410	Single 2000 x 900 x 130mm
BEM046425	Small 1920 x 820 x 130mm
BEM046460	King Single 2030 x 1050 x 130mm
BEM046475	Single - 220mm Safety Sides 2000 x 900 x 130mm

APPLICATIONS:

- Home Care
- Palliative Care
- Adolescents
- Small adults
- Users with chronic pain considerations

EFFECTIVE WEIGHT RANGE:
30-120KG

LIFECOMFORT® DELUXE



The **LIFECOMFORT® DELUXE** includes three generous foam layers; visco-elastic memory foam upper; high resilience comfort layer; and a supportive base. An ideal solution for Aged Care it provides users with an immersive, pressure relieving support surface for lasting comfort.



FOAM LAYERS

- IMMERSIVE SURFACE**
Immersive visco-elastic
- COMFORT LAYER**
High resilience foam for support and pressure redistribution
- SAFETY & SUPPORT BASE**
Strong base support for maximum durability


APPLICATIONS:

- Aged Care
- Ambulant patients
- Home Care


BEM046510	Single 2000 x 900 x 150mm
BEM046525	Small 1920 x 820 x 150mm
BEM046455	King Single 2030 x 1050 x 150mm



LIFECOMFORT® PRESSURE REDUCING



The **LIFECOMFORT® PRESSURE REDUCING** does exactly that. With visco-elastic memory foam upper, plush comfort and resilient support layers and a strong base it provides four tiers of immersive comfort and support for users requiring additional pressure relief. Integrated sidewalls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.



PROFILE HINGES



FOAM LAYERS

- IMMERSIVE SURFACE**
Immersive visco-elastic
- COMFORT LAYER**
Plush support layer catering to smaller, lighter weight users
- SUPPORT LAYER**
High resilience foam for support and pressure redistribution
- SAFETY & SUPPORT BASE**
Strong base support for maximum durability
- INTEGRATED SIDE WALLS**
Strong rigid side layers for transfer efficiency and safety

APPLICATIONS:


- Aged Care
- Ambulant patients
- Home Care
- Palliative Care

BEM046610	Single 2000 x 900 x 150mm	BEM046675	Double 2030 x 1370 x 150mm
BEM046625	Small 1920 x 820 x 150mm	BEM046680	Queen* 2030 x 1520 x 150mm
BEM046660	King Single 2030 x 1050 x 150mm		



*Maximum user weight 300kg (150kg per user).

LIFECOMFORT® PREMIER CARE



The higher weight rating of the **LIFECOMFORT® PREMIER CARE** offers greater scope for managing larger residents and respite services. Three tiers of support layers, including an immersive visco-elastic memory foam upper layer. The reinforced hinged sidewalls are designed so the mattress closely mirrors the contouring profile of an adjustable bed. Sidewall rigidity provides stability for assisted transfers or ambulant patients bed entry or exit, and accentuates the sleep surface edge.



FOAM LAYERS

- IMMERSIVE SURFACE**
Highly breathable and immersive visco-elastic
- COMFORT LAYER**
High resilience foam for diverse user weight distribution
- SAFETY & SUPPORT BASE**
Strong base support for maximum durability
- INTEGRATED SIDE WALLS**
Strong rigid side layers for transfer efficiency and safety

APPLICATIONS:

- Aged Care
- Ambulant patients
- Home Care
- Palliative Care
- Respite Services

BEM046725	Single 2000 x 900 x 150mm
BEM046730	King Single 2000 x 1050 x 150mm



SUB-ACUTE AND HOSPITAL

FOCUSED PRESSURE REDISTRIBUTION



Engineered to suit the broad patient weight requirements and pressure care needs of acute care environments. Reinforced hinged sidewalls assist transfer efficiency and accentuate the sleep surface edge. Layered sleep surface maximises pressure redistribution and immersion for patient comfort and improved sleep quality.

- 1. IMMERSIVE SURFACES**

 - Envelops vulnerable and bony prominences.
 - Evenly distributes surface contact for shear mitigation and improved skin microclimate.
 - Contours and cradles around the user for personalised pressure relief and comfort.
- 2. COMFORT LAYERS**

 - Provides an additional layer of immersive comfort and envelopment.
 - Maintains postural integrity when profiling with an adjustable bed.
- 3. SUPPORT LAYERS**

 - Allows upper layers to actively compress, maximising surface and user contact.
 - Supportive barrier for enhanced comfort in the seated/fowler position.
- 4. INTEGRATED SIDE WALLS**

 - Improves bed mobility offering a firm base of support for entry and exit.
 - Side wall firmness increases bed edge awareness for improved falls prevention.


CLINICAL BENEFITS

Shear Mitigation 	Friction Reduction 	Skin Microclimate
Entry & Exit Safety 	Immersion 	Higher SWL

LIFECOMFORT® ACUTE CARE



The **LIFECOMFORT® ACUTE CARE** mattress has three foam layers including a castellated patient surface providing users with immersive, pressure relieving support and air flow for enhanced skin microclimate. The U-shaped base and sidewall structure provides support to assist with patient transfers as well as creates roll resistance at the sleep surface edges.



FOAM LAYERS

- IMMERSIVE CASTELLATED SURFACE**
Castellated cells offer greater airflow and transitional support
- COMFORT LAYER**
Visco-elastic comfort and envelopment
- SAFETY & SUPPORT BASE**
Strong U-shaped base for maximum durability and assistive support
- INTEGRATED SIDE WALLS**
Strong rigid side layers for transfer efficiency and safety

APPLICATIONS:

- Hospital
- Aged Care
- Home Care
- Palliative Care


BEM046615	Single 1980 x 880 x 150mm
BEM046616	King Single 2030 x 1050 x 150mm



LIFECOMFORT® BARIATRIC



The **LIFECOMFORT® BARIATRIC** is a robust support surface designed for larger users up to 300kg. Four foam tiers with a visco-elastic memory foam upper provide immersion for pressure reduction and comfort as well as support. Side wall foam density is strengthened to assist heavier users for bed entry and exit.



PROFILE HINGES



FOAM LAYERS

- IMMERSIVE SURFACE**
Immersive visco-elastic
- COMFORT LAYER**
Transition layer that boosts visco layer efficiency under higher weight loads
- SUPPORT LAYER**
High resilience foam for diverse user weight and distribution
- SAFETY & SUPPORT BASE**
Strong base support for maximum durability
- INTEGRATED SIDE WALLS**
Strong rigid side layers for transfer efficiency and safety

APPLICATIONS:

- Hospital
- Aged Care
- Home Care
- Palliative Care
- Bariatric Care


BEM046700	King Single 2030 x 1050 x 200mm
-----------	------------------------------------



LIFECOMFORT® MAXI LIFESTYLE



The **LIFECOMFORT® MAXI LIFESTYLE** is a robust support surface designed for larger users up to 350kg. Three foam tiers provide immersion for pressure reduction and comfort as well as support. Side wall foam density is strengthened to assist heavier users for bed entry and exit. Fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care, Aged Care and up to Acute Care environments.



PROFILE HINGES



FOAM LAYERS

- COMFORT LAYER**
Transition layer that boosts visco layer efficiency under higher weight loads
- SUPPORT LAYER**
High resilience foam for diverse user weight and distribution
- SAFETY & SUPPORT BASE**
Strong base support for maximum durability
- INTEGRATED SIDE WALLS**
Strong rigid side layers for transfer efficiency and safety

APPLICATIONS:

- Hospital
- Aged Care
- Home Care
- Palliative Care
- Bariatric Care

BEM046710	Super King 2060 x 1150 x 200mm
-----------	-----------------------------------

EFFECTIVE
WEIGHT RANGE:
150-350KG

LIFECOMFORT® PREMIER CARE ACUTE



The **LIFECOMFORT® PREMIER ACUTE CARE** has a state of the art immersive memory foam layer up to 8x more breathable than conventional memory foam for an optimal sleep temperature and skin microclimate management. Greater foam elasticity makes repositioning and moving on and off the mattress easier. The multi-layered support surface with castellated cells, closely contours to the user's body for optimal pressure redistribution and comfort.



PROFILE HINGES



FOAM LAYERS

- IMMERSIVE SURFACE**
Ultra breathable Memory Foam for superior immersion
- CASTELLATED DUAL COMFORT LAYER**
Castellated cells for greater airflow and transitional support
- SAFETY & SUPPORT BASE**
Strong base support for maximum durability
- INTEGRATED SIDE WALLS**
Strong rigid side layers for transfer efficiency and safety

APPLICATIONS:

- Hospital
- Aged Care
- Home Care
- Palliative Care
- Respite Services

BEM046735	Single 2000 x 900 x 150mm
BEM046740	King Single 2000 x 1050 x 150mm


EFFECTIVE
WEIGHT RANGE:
30-250KG

*Up to 8x more breathable than conventional memory foams
aspire LIFECOMFORT® Foam Mattresses

HYBRID COMFORT



The **HYBRID COMFORT** is an Acute mattress system offering foam-in-cell technology enabling flexible operation as a high comfort foam static support surface or a powered-up alternating pressure support surface. Specifically re-engineered for increased patient comfort, foam layering now effectively supports patient weight ranging from 50-230kg.



FEATURES AND BENEFITS

Immersive Visco-elastic Memory Foam Surface: Top layer is now twice as thick to offer significant improvement in user comfort

Intelligent Hose Management: When not in use the air-hose can be stored internally, reducing infection control and OH&S risks.

Integrated Shear Reduction: Memory foam layer moves independently of the cell surface reducing shear when bed is adjusted.

Optimised Heel Protection: Tapered heel zone redistributes pressure to lower leg, away from the vulnerable heel region. The castellated surface collapses in a wave-like movement reducing shear.

Innovative pump technology: Reduces noise and vibration for uninterrupted sleep.

APPLICATIONS:

- Community
- Aged Care
- Hospital
- Ambulant patients
- Non-ambulant patients

PTM638005	Single 1980 x 895 x 170mm
PTM638006	King Single 1980 x 1050 x 170mm



Pump
PTM638030



EFFECTIVE
WEIGHT RANGE:
50-230KG

RELATED PRODUCTS



MATTRESS EXTENSION BOLSTERS

- Addresses variability of bed lengths reducing the risk of injury to the patient.

BEM046762	Single - Short - 900 900 x 100 x 130mm
BEM046764	Single - Short - 1050 1050 x 100 x 130mm
BEM046772	King Single - Long - 900 900 x 200 x 130mm
BEM046774	King Single - Long - 1050 1050 x 200 x 130mm



LIFECOMFORT FALL SAFETY MAT

- Prevents serious injury to patient in the event of a fall from bed.
- Nylon cover with reflective strip for easier visibility in low light settings.

BEA008550	Fall Safety Mat 1760 x 600 x 50mm
-----------	--------------------------------------



LIFECOMFORT EVACUATION SHEET

- Designed for emergency evacuation of non-ambulant patients.
- Fits discreetly underneath the resident's mattress and remains readily accessible.
- Fire retardant and highly durable.

TSE672526	Evacuation Sheet 1900 x 900mm (S.W.L. 300kg)
-----------	---



LIFECOMFORT PRESSURE SAFETY UNDERLAYS

- Provides "bottoming out" protection for air alternating overlays.
- Vapour permeable (breathable), antibacterial and water resistant, with a fire retardant outer cover.

BEM047050	Single 2000 x 900 x 50mm (S.W.L. 135kg)
BEM047060	King Single 2030 x 1050 x 50mm (S.W.L. 135kg)



BED RAIL PROTECTORS

- Protects patients from injury by preventing limbs from protruding through bed rails as well as any contact.
- Velcro fastening straps allows for easy application and removal.
- Easy to clean Polyurethane fabric to assist with infection control management.

BEA010775	Small (Pair) 1260 x 670mm
BEA010776	Large (Pair) 1580 x 980mm



07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034

sales@powermobility.com.au

powermobility.com.au



ABBRT1.43 1110211

1300 133 120
www.aspirecare.com.au

 **aspire LIFEComfort®**