



# WHILL C2 CUSTOM BACKREST CANE ASSEMBLY

## INSTALLATION

Installing these canes is quick and easy, with adjustable width of 16-20" provided for aftermarket backrests to be fitted.

The standard configuration stands at 300mm in height, accommodating most common backrest heights. If there is a need for a taller backrest, this design allows for a seamless upgrade: simply replace the vertical upright with a 500mm tube to increase the height of the backrest cane assembly.

## PRODUCT FEATURES

- Suits most aftermarket backrests
- Width adjustable - Suits 16 - 20" backrests
- Height adjustable – Taller backrest

upright allows installation of taller custom backrest cushions

**Detailed Fitment  
Instructions Video**  
for the custom backrest  
cane assembly



## BACKREST

Code	Description
WCA870456	Custom Backrest Frame - Standard
WCA870457	Custom Backrest Frame - Tall

## SUITS CHAIRS

Code	Description
PWS655400	Whill C2 - 450mm (18") Power Wheelchair
PWS655410	Whill C2 - 500mm (20") Power Wheelchair