



DASH WHEELCHAIR

The Aspire Dash is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends.

Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support. The 135kg safe working load caters to a broad range of users.



1 YEAR WARRANTY
Components

2 YEAR WARRANTY
Frame



UPHOLSTERY
Orange/Black

FRAME
Silver



UPHOLSTERY
Green/Black

FRAME
Silver

FEATURES & BENEFITS

- Available as attendant propelled or self-propelled
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers
- Contoured and padded seating improves comfort and postural support
- Attendant operated hand brakes with parking lock
- Simple and intuitive folding for easy transportation
- Removable leg rests and folding back for compact storage
- Rear storage pocket for essentials
- Anti-tippers available (WCA848270)

Self-Propelled



Attendant Propelled



OVERALL LENGTH	OVERALL WIDTH	OVERALL HEIGHT	SEAT DEPTH	SEAT WIDTH	ANTERIOR SEAT HEIGHT	POSTERIOR SEAT HEIGHT	BACKREST HEIGHT	ARMREST HEIGHT	REAR WHEEL	MIN. CHAIR WEIGHT	TOTAL CHAIR WEIGHT
(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(inch)	(inch)	(kg)	(kg)

MWS449850 Attendant Propelled Orange	995	655	930	420	450	495	445	430	200	13"	8.6	9.6
MWS449855 Attendant Propelled Green	995	655	930	420	450	495	445	430	200	13"	8.6	9.6
MWS449860 Self Propelled Orange	1030	685	930	420	450	495	445	430	200	22"	10.7	11.7
MWS449865 Self Propelled Green	1030	685	930	420	450	495	445	430	200	22"	10.7	11.7

ASFY77.3 020322

1300 133 120
www.aspirecare.com.au

power: mobility
The power to move you

07 3265 4663
2/452 Bilsen Rd, Geebung QLD 4034
sales@powermobility.com.au • powermobility.com.au