



# SOCIALITE WHEELCHAIR

*Arrive in style*

The Aspire Socialite is an ultra light weight transport chair perfect for short trips and outings with friends. The generous 115kg safe working load caters to a broad range of users.

Padded and breathable seat and back provide optimum comfort and postural support. Cut away arms improve table access and the offset uprights create hip space for coats and warmer clothing.

A photograph showing an elderly man with white hair sitting in a red and black transport wheelchair. He is wearing a light blue and white checkered shirt and blue jeans. A woman in a dark jacket is pushing the wheelchair from behind. In the background, another woman in a blue patterned shirt is standing and talking to him. The setting appears to be a kitchen or a dining area with a white tiled wall and a stainless steel stove.

LESS THAN  
**9KG**

# FEATURES & BENEFITS

- Sleek European design
- Stylish contemporary colours
- One of the lightest transport chairs available
- Perfect for outings and appointments
- Attendant operated hand brakes with parking lock
- Simple and intuitive folding for easy transportation
- Removable leg rests and folding back for compact storage
- Contoured seating improves comfort and postural support
- Rear storage pocket for essentials



UPHOLSTERY  
Red/Black

FRAME  
Red



UPHOLSTERY  
Silver/Black

FRAME  
Silver

**1** YEAR WARRANTY  
Components

**2** YEAR WARRANTY  
Frame



| ALUMINIUM   | OVERALL WIDTH (mm) | SEAT DEPTH (mm) | SEAT WIDTH (mm) | SEAT HEIGHT ANTERIOR (mm) | SEAT HEIGHT POSTERIOR (mm) | BACKREST HEIGHT (mm) | ARMREST HEIGHT (mm) | PRODUCT WEIGHT (kg) | MIN. CHAIR WEIGHT (kg) |
|---|--------------------|-----------------|-----------------|---------------------------|----------------------------|----------------------|---------------------|---------------------|------------------------|
| <b>Socialite</b><br>MWS449800 - Silver<br>MWS449810 - Red | 595                | 405             | 450             | 475-500                   | 440                        | 410                  | 200                 | 8.9                 | 8.1                    |



1300 133 120  
www.aspirecare.com.au



07 3265 4663

2/452 Bilsen Rd, Geebung QLD 4034  
sales@powermobility.com.au • powermobility.com.au