



VIDA WHEELCHAIR

The Aspire VIDA is an extremely lightweight wheelchair that can be lifted in and out of a vehicle easily by carers. The VIDA range comes in various fun, bold colours and is ideal for quick trips to the shops, appointments, and catching up with friends. The ergonomic backrest and seat shape mimic body contours for comfort and postural support. The high safe working load (135kg / 21st) caters to a wide range of users.



UPHOLSTERY
Mustard/Black
FRAME
Grey
SEAT WIDTH
400mm



UPHOLSTERY
Teal/Black
FRAME
Grey
SEAT WIDTH
450mm



UPHOLSTERY
Magenta/Black
FRAME
Grey
SEAT WIDTH
450mm



UPHOLSTERY
Olive/Black
FRAME
Grey
SEAT WIDTH
500mm

Minimum
Chair Weight
from
7.5kg



5 YEARS WARRANTY
Frame

1 YEAR WARRANTY
Components

FEATURES & BENEFITS

- Available as attendant propelled or self-propelled
- Three seat width sizes: 400mm, 450mm, and 500mm, each with its own unique colour
- The 400mm seat width comes in Mustard
- The 450mm seat width is available in Teal and Magenta
- The 500mm seat width is available in Olive
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers
- Contoured and padded seating improves comfort and postural support
- Attendant operated hand brakes with parking lock
- Simple and intuitive folding for easy transportation
- Removable leg rests and folding back for compact storage
- Rear storage pocket for essentials
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)



Attendant Propelled (AP) & Self Propelled (SP) Vida Specifications

400mm Mustard



MWS449874 - AP
MWS449876 - SP

450mm Teal or Magenta



MWS449878 - AP
MWS449880 - SP

500mm Olive



MWS449886 - AP
MWS449888 - SP



VIDA	Overall Length mm/in	Overall Width mm/in	Overall Height mm/in	Seat Depth mm/in	Seat Width mm/in	Anterior Seat Height mm/in	Posterior Seat Height mm/in	Backrest Height mm/in	Armrest Height mm/in	Min Chair Weight kg/lbs	Rear Wheel in
Attendant Propelled											
MWS449874	995/39	600/23.5	930/36.5	420/16.5	400/16	500/19.5	460/18	430/17	200/8	8.8/19.5	13
MWS449878	995/39	650/25.5	930/36.5	420/16.5	450/18	500/19.5	460/18	430/17	200/8	8.8/19.5	13
MWS449882	995/39	650/25.5	930/36.5	420/16.5	450/18	500/19.5	460/18	430/17	200/8	8.8/19.5	13
MWS449886	995/39	680/27	930/36.5	420/16.5	500/20	500/19.5	460/18	430/17	200/8	9.3/20.5	13
Self Propelled											
MWS449876	1040/41	630/25	950/37.4	420/16.5	400/16	500/19.5	460/18	430/17	200/8	7.5/16.5	22
MWS449880	1040/41	690/27	950/37.4	420/16.5	450/18	500/19.5	460/18	430/17	200/8	7.5/16.5	22
MWS449884	1040/41	690/27	950/37.4	420/16.5	450/18	500/19.5	460/18	430/17	200/8	7.5/16.5	22
MWS449888	1040/41	710/28	950/37.4	420/16.5	500/20	500/19.5	460/18	430/17	200/8	8/17.6	22

*Range of motion dimensions are accurate to a +/-10mm tolerance