

	Scooter Code	Image
COMPACT	<b>BOOT SCOOTER</b> (SCT670613)	
	<b>GO-GO® LX + CTS SCOOTER</b> (SCT670612)	
	<b>SONIC SCOOTER</b> (SCT670665)	
	<b>QUEST FOLDING COMPACT SCOOTER</b> (SCT670630)	
	<b>BOOT SCOOTER SUPALITE</b> (SCT670616)	
MEDIUM	<b>GT SCOOTER</b> (SCT670670)	
	<b>PATHRIDER® 10 SCOOTER</b> (SCT670615)	
	<b>PATHRIDER® 10 DX SCOOTER</b> (SCT670620)	
	<b>VICTORY 10 LX SCOOTER</b> (SCT670640)	
	<b>GOPHER DX SPORT SCOOTER</b> (SCT670638)	
LARGE	<b>PATHRIDER® 130XL PURSUIT SCOOTER</b> (SCT670571)	
	<b>PATHRIDER® 140 XL SCOOTER</b> (SCT670576)	
	<b>CELEBRITY® DX XL SCOOTER</b> (SCT670540)	
	<b>SPORTRIDER SCOOTER</b> (SCT670585)	
	<b>RANGERIDER SCOOTER</b> (SCT670581)	

Terrain	Battery	Max Range Per Charge*	Safe Working Load
Level & smooth terrain. Short distances.	2x 12V 12AH	11.5 km	118 kg
Level & smooth terrain. Short distances.	2x 12V 18AH	17 km	136 kg
Level & smooth terrain. Short distances.	2 x 12V 20AH	20 km	136 kg
Able for more varied terrain. Short to medium distances.	1x 24V 11.5AH Li-Polymer	20 km	115 kg
Level & smooth terrain. Short distances.	1x 24V 12AH	14.4 km	120 kg
Able for more varied terrain. Short to medium distances.	2 x 12V 35AH	30 km	159 kg
Able for more varied terrain. Short to medium distances.	2x 12V 31AH (40AH optional)	40 km	159 kg
Able for more varied terrain. Short to medium distances.	2x 12V 31AH (40AH optional)	33 km	159 kg
Can handle rougher terrain well. Long distances.	2x 12V 40AH	26 km	181 kg
Able for more varied terrain. Short to medium distances.	2x 12V 50AH	33 km	160 kg
Can handle rougher terrain well. Long distances.	2x 12V 55AH	25 km	159 kg
Can handle rougher terrain well. Long distances.	2x 12V 70AH	35 km	181 kg
Can handle rougher terrain well. Long distances.	2x 12V 40AH	40 km	159 kg
Can handle rougher terrain well. Long distances.	2x 12V 70AH	50 km	181 kg
Can handle rougher terrain well. Long distances.	2x 100AH	30 km	159 kh

Features
- Up to 6km/h - Disassembles for easy transport
- Up to 6.4km/h - Disassembles for easy transport - Comfort Trac front & rear suspension (CTS)
- Up to 6.4km/h - Disassembles in easy steps
- Up to 8.8km/h - Great outdoor stability, ease of use and maximum value
- Up to 6.5km/h - Disassembles for easy transport - Easy-to-grip tiller adjustment knob
- Up to 10km/h - Front & rear suspension
- Up to 9km/h - Comfort Trac front & rear suspension (CTS)
- Up to 8.5km/h - Delivers high performance and feather-touch disassembly - Can upgrade to 40AH batteries (x2) for improved performance
- Up to 9km/h - Front & rear suspension - Can upgrade to 70AH batteries (x2) for improved performance
- Up to 10km/h - Front & rear suspension - High visibility, front and rear lighting package
- Up to 10km/h - Front & rear suspension - Can upgrade to 100AH batteries (x2) for improved performance
- Up to 10km/h - Front & rear suspension - Can upgrade to 100AH batteries (x2) for improved performance
- Up to 10km/h - Front & rear suspension - Can upgrade to 100AH batteries (x2) for improved performance
- Up to 10km/h - Front & rear suspension - CTS Suspension includes adjustable shocks for greater comfort - Backup sensors with LED light

\*Range Per Charge varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition & tyre condition. Battery upgrades (where possible) improve Range Per Charge.

## SCOOTER MOBILITY

At the heart of Pride's philosophy is the understanding that mobility products empower lives. Through industry-leading, consumer-inspired designs, Pride offers an innovative, distinctive and performance driven comprehensive line of mobility products to help our customers lead the most active, independent lifestyles.

Aidacare works in close association with Pride Australia to promote and support Pride mobility products

- Pride Scooter products are all registered in the ARTG
- Please always read and understand your owners handbook before operating your new mobility scooter
- Please take care to observe safety and the speed limit requirements of your state when driving your mobility scooter
- Specification and other details on this information sheet may be subject to change without notice

- 1 Comfortable front-mounted post seat
- 2 Back-lit battery gauge on monitor
- 3 Front basket as standard
- 4 Sleek & sporty design
- 5 LED headlight, tail light and indicators on some models
- 6 Low profile, puncture proof tyres or pneumatic tyres



### SAFE USE VIDEO

At Aidacare, we value the safety of our customers above all else. There are a few safety considerations we'd like you to understand before operating your device.

We strongly recommend all new mobility device customers watch our 7 minute YouTube video, 'Getting to Know your Power Mobility Device' that covers the above 8 points in detail and more.



## MOBILITY SCOOTER SELECTION GUIDE



# MOBILITY SCOOTER SELECTION GUIDE

## aspire BOOT SCOOTER SCT670613

### FEATURES

- Feather touch disassembly into 5 pieces for easy transportation & storage
- Comfortable front-mounted post seat with arm rests for stability
- Easy to remove rechargeable battery
- Soft and easy to adjust tiller
- Modular design for easy servicing

### SPECIFICATIONS

Length: 1003 mm      Colour: Red  
 Width: 495 mm  
 Product Weight: 36.6 kg  
 Maximum Speed: 6.4 km/h  
 Range: 11.5 km    HP: 1.6  
 Battery Size: 2 x 12V (12AH)



118kg  
S.W.L

## QUEST FOLDING COMPACT SCOOTER SCT670630

### FEATURES

- Sleek, stylish and extremely easy to operate
- Ready for transport and folded in just 3 easy steps
- Can be stowed upright or laid down
- Great performance both indoors and out
- Articulated front end
- Front and rear LED headlights

### SPECIFICATIONS

Length: 930 mm      Colour: As Shown  
 Width: 485 mm  
 Product Weight: 23.5 kg  
 Maximum Speed: 6 km/h  
 Range: 20 km    HP: 0.362  
 Battery Size: 1 x 24V (11.5 AH) Li-Polymer



115kg  
S.W.L

## PATHRIDER® 10 DX SCOOTER SCT670620

### FEATURES

- High performance sleek and sporty mobility scooter
- Low profile, puncture proof tyres
- Comfortable front-mounted post seat with armrests for stability
- Back-lit battery gauge on monitor
- Rugged front and rear suspension

### SPECIFICATIONS

Length: 1193 mm      Colour: Blue  
 Width: 565 mm  
 Product Weight: 85 kg  
 Maximum Speed: 9 km/h  
 Range: 21 / 25 km    HP: 1.8  
 Battery Size: 2 x 12V (31AH) (40AH optional)



159kg  
S.W.L

## PATHRIDER® 140 XL SCOOTER SCT670576

### FEATURES

- Large 13" pneumatic tyres
- Luxurious suspension seat with arm rests, head rest and lumbar support
- Front and rear suspension and sealed hydraulic brakes
- Automatic full lighting system
- Back lift battery gauge
- Wrap around easy drive tiller

### SPECIFICATIONS

Length: 1428 mm      Colour: Red  
 Width: 692 mm      Colour: Yellow  
 Product Weight: 118 kg  
 Maximum Speed: 10 km/h  
 Range: 35 km    HP: 4.6  
 Battery Size: 2 x 12V (70AH)



181kg  
S.W.L

## GO-GO® LX + CTS SCOOTER SCT670612

### FEATURES

- CTS Suspension (Comfort-Trac Suspension) brings a new level of performance to travel mobility
- Lightweight mobility scooter
- Disassembles into 5 pieces for easy transportation and storage
- Easy to remove rechargeable battery
- Soft and easy to adjust tiller

### SPECIFICATIONS

Length: 1040 mm      Colour: Red  
 Width: 520 mm  
 Product Weight: 39.5 kg  
 Maximum Speed: 6.4 km/h  
 Range: 17 km    HP: 1.6  
 Battery Size: 2 x 12V (18AH)



136kg  
S.W.L

## aspire GT SCOOTER SCT670670

### FEATURES

- Offers full advanced suspension for a smoother, more comfortable ride over varied terrain
- Full advanced suspension
- Front and rear LED lighting illuminates your path
- Deluxe tiller with USB charger and cup holder

### SPECIFICATIONS

Length: 1200 mm      Colour: Red  
 Width: 605 mm  
 Product Weight: 82 kg  
 Maximum Speed: 10 km/h  
 Range: 30 km    HP: 0.536  
 Battery Size: 2 x 12V (35AH)



159kg  
S.W.L

## VICTORY 10 LX SCOOTER SCT670640

### FEATURES

- CTS Suspension (Comfort-Trac front and rear suspension)
- Sleek, modern design
- Feather-touch disassembly permits simple frame separation with only one hand
- Full LED lighting package that includes curb and rear lights and directional signals

### SPECIFICATIONS

Length: 1194 mm      Colour: Red  
 Width: 565 mm  
 Product Weight: 71 kg  
 Maximum Speed: 10 km/h  
 Range: 26 km    HP: 0.570  
 Battery Size: 2 x 12V (40AH)



181kg  
S.W.L

## CELEBRITY® DX XL SCOOTER SCT670540

### FEATURES

- Large knobby tyres
- Deluxe front-mounted high back swivel seat with headrest and arm rests for stability
- Lumbar support and manual recline
- State of the art rear suspension
- Unique easily adjusted steering mount with modern LED lighting

### SPECIFICATIONS

Length: 1300 mm      Colour: Red  
 Width: 640 mm  
 Product Weight: 91 kg  
 Maximum Speed: 10 km/h  
 Range: 40 km    HP: 2.83  
 Battery Size: 2 x 12V (40AH)



159kg  
S.W.L

## aspire BOOT SCOOTER SUPALITE SCT670616

### FEATURES

- Compact, foldable seat with armrests
- Easy-to-grip tiller adjustment knob
- Stable, front frame-mounted seat post
- A standard front basket provides storage
- Easy to transport at home, hotels, or airports

### SPECIFICATIONS

Length: 1022 mm      Colour: White  
 Width: 490 mm  
 Product Weight: 8.5 kg  
 Maximum Speed: 6.5 km/h  
 Range: 14.4 km    HP: 8.9  
 Battery Size: 24V (12AH)



120kg  
S.W.L

## GOPHER DX SPORT SCT670638

### FEATURES

- Speed reduction technology automatically reduces speed when turning corners
- Independent front suspension, mono shock rear suspension
- Easy-drive tiller with wraparound handles
- High visibility, front and rear lighting package

### SPECIFICATIONS

Length: 1370 mm      Colour: Blue  
 Width: 670 mm  
 Product Weight: 102.5 kg  
 Maximum Speed: 10 km/h  
 Range: 33 km    HP: 20.62  
 Battery Size: 12V (50AH)



160kg  
S.W.L

## aspire SONIC SCOOTER SCT670665

### FEATURES

- Quickly dismantles into five light weight pieces for easy transport and storage, no tools required
- Portable and durable scooter
- Puncture-proof tyres
- Detachable battery pack
- Front and rear LED light
- Easy access battery compartment
- Adjustable angle tiller with delta handle

### SPECIFICATIONS

Length: 1040 mm      Colour: Red  
 Width: 495 mm  
 Product Weight: 45 kg  
 Maximum Speed: 6.8 km/h  
 Range: 20 km    HP: 0.241  
 Battery Size: 2 x 12V 20AH



136kg  
S.W.L

## PATHRIDER® 10 SCOOTER SCT670615

### FEATURES

- High performance, sleek and sporty mobility scooter
- Low profile, puncture proof tyres
- Back lift battery gauge
- Comfortable front-mounted post seat with armrests for stability
- Eight-position seat rotation
- Front basket as standard
- Wrap around easy drive tiller

### SPECIFICATIONS

Length: 1190 mm      Colour: Blue  
 Width: 565 mm  
 Product Weight: 63.5 kg  
 Maximum Speed: 8.5 km/h  
 Range: 21 / 25 km    HP: 1.8  
 Battery Size: 2 x 12V (31AH) (40AH optional)



159kg  
S.W.L

## aspire PATHRIDER® 130XL PURSUIT SCOOTER SCT670571

### FEATURES

- High performance, sleek and sporty
- Large 13" pneumatic tyres
- Luxurious suspension seat with arm rests, head rest and lumbar support
- Rugged front and rear suspension
- Automatic full lighting system
- Long lasting high intensity LED kerb light, headlight, tail light and indicators

### SPECIFICATIONS

Length: 1359 mm      Colour: Red  
 Width: 635 mm  
 Product Weight: 73 kg  
 Maximum Speed: 12.4 km/h  
 Range: 25 km    HP: 3.1  
 Battery Size: 2 x 12V (55AH)



159kg  
S.W.L

## SPORTRIDER SCOOTER SCT670585

### FEATURES

- A sporty, sleek outdoor scooter
- Full digital displaying speedometer/odometer
- Full suspension system
- High back, adjustable seat
- Powerful front/rear lighting kit
- Large rear storage pod included as standard

### SPECIFICATIONS

Length: 1650 mm      Colour: Black  
 Width: 780 mm  
 Product Weight: 120.5 kg  
 Maximum Speed: 12 km/h  
 Range: 50 km    HP: 4.1  
 Battery Size: 2 x 12V (70AH)



181kg  
S.W.L

## RANGERIDER SCT670581

### FEATURES

- Limited recline, black vinyl, high-back seat with lumbar and pillow-top headrest, flip-up height, and angle adjustable armrests
- CTS Suspension includes adjustable shocks for greater comfort
- Backup sensors with LED light
- Dual hydraulic brakes

### SPECIFICATIONS

Length: 1560 mm      Colour: Stone Grey  
 Width: 830 mm  
 Maximum Speed: 10 km/h  
 Range: 30 km    HP: 18.6  
 Battery Size: (2) (100AH)



159kg  
S.W.L