



# VOGUE WALKER RANGE

aspire for...

Comfort

Stability

Style



[www.aspirecare.com.au](http://www.aspirecare.com.au)

1300 133 120



# VOGUE CARBON FIBRE SEAT WALKER



One of the lightest seat walkers on the market! Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super lightweight carbon fibre frame for a superior walking experience.



## MODEL

WAF705350 | Vogue Carbon Fibre - Medium - Black  
 WAF705360 | Vogue Carbon Fibre - Tall - Candy Red

## FEATURES & BENEFITS shared by the Vogue Carbon Fibre and Lightweight Aluminium Seat Walkers



### European Design

Ergonomic design offers personalised height adjustment, promoting correct posture



### Compact

Turning circle for improved mobility in tight spaces



### Simple

Folding mechanism locks for easy transportation and storage



### Handbrakes

Lock for user safety and stability

	LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT** (kg)	BAG S.W.L (kg)
MED	680	600	695 - 920	240	460	550	225	8"	4.9	6.2	5
TALL	680	605	770 - 1040	240	460	620	225	8"	5.2	6.3	5

# VOGUE LIGHTWEIGHT SEAT WALKER



This modern seat walker has a clean and crisp frame design, making it a standout when going for a walk or browsing the shops!



Luxury Look & Feel  
Danish Designed

AVAILABLE COLOURS:



Red



Champagne

WEIGHS  
JUST OVER  
5KG



Discreet Cable Management



**Lightweight Aluminium Body**  
Improves manoeuvrability and portability

## MODEL

WAF705310RE | Vogue Lightweight Aluminium - Red  
WAF705310CH | Vogue Lightweight Aluminium - Champagne



### Storage

Removable bag carries up to 5kg of essentials



### Backrest

Deep seat and wide backrest for improved comfort and support



### Bonus

Walking cane holder

**1 YEAR LIMITED WARRANTY**

AS/NZ Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

OVERALL			SEAT		FOLDED	WHEEL	MINIMUM FUNCTIONAL	TYPICAL FUNCTIONAL	BAG
LENGTH (mm)	WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	WIDTH (mm)	HEIGHT (mm)	WIDTH (mm)	WEIGHT (kg)	WEIGHT** (kg)	S.W.L (kg)
655	620	785 - 920	240	465	540	220	5.3	6.4	5

\* Carbon Fibre Walker - Medium  
\*\* Includes accessories, bag and backstrap.



Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super sized front wheels can manage any surface, perfect for the adventurous user!

**Wider Seat**  
supports up to  
200kg

**Folds to**  
245mm

**Large**  
Multi-Terrain  
Wheels



Luxury Look & Feel  
Danish Designed

AVAILABLE COLOURS:



Black

**AutoLocate**

Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate.



**MODEL**

WAF705450 | Vogue Adventure - Black

**FEATURES & BENEFITS**



**European Design**

Ergonomic design offer personalised height adjustment, promoting correct posture



**Compact**

Turning circle for improved mobility in tight spaces



**Simple**

Folding mechanism locks for easy transportation and storage

**1 YEAR LIMITED WARRANTY**

✓ AS/NZ Approved

✓ Vogue is registered in the ARTG, AUST R ARTG 271255



**Handbrakes**

Lock for user safety and stability



**Storage**

Removable bag carries up to 5kg of essentials



**Bonus**

Walking cane holder

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	SEAT HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)
770	695	900 - 1035	230	540	530	245	Front 11" Rear 9"	8.8	10	5



For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.



Luxury Look & Feel  
Danish Designed

AVAILABLE COLOURS:



White / Black



**Lightweight Aluminium Body**  
Improves portability and manoeuvrability



MODEL

WAF705400 | Vogue Forearm - White/Black

FEATURES & BENEFITS



**European Design**  
Ergonomic design offer personalised height adjustment, promoting correct posture



**Simple**  
Folding mechanism locks for easy transportation and storage



**Bonus**  
Walking cane holder

**1 YEAR LIMITED WARRANTY**

✓ AS/NZ Approved

✓ Vogue is registered in the ARTG, AUST R. ARTG 271255



**Handbrakes**  
Lock for user safety and stability



**Storage**  
Removable bag carries up to 5kg of essentials

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)
660	620	1020 - 1270	240	460	620	225	8"	9.1	10*	5

\* Includes accessories, bag and backstrap.



Ideal for use within the home, the Indoor Walker is perfect for portable dining and in home mobility assistance.




Luxury Look & Feel  
Danish Designed

AVAILABLE COLOURS:



Red / Black



Champagne / Black



Folds to  
270mm



Portable Dining

Safely collect and transport meals or hot drinks

On-cart Dining



Turn Around.  
Slide Over Knee.  
Lock & Enjoy.



MODEL

- WAF705100RE | Vogue Indoor - Red/Black
- WAF705100CH | Vogue Indoor - Champagne/Black

FEATURES & BENEFITS



**European Design**  
Ergonomic design offer personalised height adjustment, promoting correct posture



**Compact**  
Turning circle for improved mobility in tight spaces



**Conveniently**  
Stands up when folded for flat storage



**Handbrakes**  
Lock for user safety and stability



**Storage**  
Removable bag carries up to 5kg of essentials



**Bonus**  
Walking cane holder

**1 YEAR LIMITED WARRANTY**

AS/NZ Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)	TRAY S.W.L (kg)
650	540	860 - 970	270	6"	4.4	5.1	5	5

\* Includes accessories, bag and tray.



Oxygen  
Bottle Holder Bag  
WAA695600



Night Light - LED Lamp  
WAA695630



Single Hand Brake System\* -  
Compatible with Aspire Vogue Carbon Walkers  
WAA695640

*\*Single Hand Brake System are not compatible with the Vogue Indoor Walker*



CRUTCHES, FRAMES & STICKS



SCOOTERS



OVERBED TABLES



BED STICKS & POLES



RAMPS



AIR MATTRESSES & FOAM MATTRESSES



COMMODOES, URINALS & BED PANS



OVER TOILET AIDS



SHOWER STOOLS & CHAIRS



RAISED TOILET SEATS



SHOWER COMMODOES



BATH BOARDS & BENCHES



ADJUSTABLE BEDS



POWER WHEELCHAIRS



PRESSURE CUSHIONS



LIFT & RECLINE CHAIRS



**power mobility**  
The power to move you

**07 3265 4663**

2/452 Bilsen Rd, Geebung QLD 4034

sales@powermobility.com.au • powermobility.com.au

