

MovingStar MS 101 SF R

Folding Transportable Power Chair

If you're looking for a lightweight power chair that offers convenience, transportability and great value, the MovingStar MS 101 SF R could be just what you need!

- **Total Weight 21.5kg**
- **Easily splits into 2 parts for transportation & storage (12.2kg + 8.8kg)**
- **Patented footrest design**
- **Optional Attendant Control**
- **Folds in seconds and fits easily in a small car**
- **Ideal for around the home, at the shops, park or family outings.**



Call us today to arrange a trial or visit our showroom to talk to our experienced team.

2/452 Bilsen Rd, Geebung, 4034 | 07 3265 4663

power
mobility
The power to move you

SPECIFICATIONS – MS 101 SF R

Total Weight: 21.5 kg

Easily splits into 2 parts: (Drive Unit - 12.2kg, Frame - 8.8kg)

Lithium Battery: (1 x 24v 10.5AH 250W) 1.2kg

Colours: Black or Blue

Maximum Speed: 6km/hr

Approx. Drive Range: 15km

User Weight Capacity: 115kg

Solid Tyres: 8" Front, 12" Rear

Seat Depth: 42cm

Seat Width: 49cm

Backrest Height: 42cm

Floor to Seat Height: 50cm

Distance between armrests: 49.5cm

Overall Length: 81cm

Overall Width: 63cm

Overall Height: 91.5cm

Folded Size: 93cmL x 63cmW x 28cmH (split)

Motor Brushless: DC 200W x 2pcs