



The innovative wheelchair wheels from trivida®

trivida® Introducing a new state-of-the-art wheelchair wheel

A triumph for freedom and independence

trivida® is the world's first wheel that makes it possible to transfer from and into the wheelchair without any obstacles. Thanks to the revolutionary triatec® technology, the upper segment of the wheel can be easily removed so that a lateral transfer from the wheelchair to the external seat and vice versa is possible without any obstacles.

Simple assembly and maintenance Get personal advice

The trivida® wheel has a universal quick-change axle and is suitable for almost all types of wheelchairs. The standard 50 mm long bush with 12 mm diameter normally guarantees installation in seconds. For wheelchairs with 1/2" axle system the wheel can also be easily installed using the bushes supplied. For the few special cases with different axle holders or wheel fastening systems we will find a customised technical solution for you. We are happy to advise you. We also recommend an annual maintenance service.

Find out more at:
www.trivida-info.com

Sales & Service Europe:
P + L Innovations GmbH
Am Krozinger Weg 11
79189 Bad Krozingen
Germany

Tel. +49 (0) 7633 920 3939
sales@trivida-info.com



The universal quick-change axle fits most types of wheelchair.



trivida®
Barrier-free living.

a P+L Innovations GmbH brand



trivida®
plastic wheel



trivida®
carbon wheel

German engineering leads the world yet again. This time, we've reinvented the wheel.

With its revolutionary triatec technology, the trivida® wheel is manufactured in Germany to the world's highest standards using high-quality, durable materials. The modern, functional design of this wheelchair wheel is groundbreakingly brilliant. It's such a deceptively simple idea, yet so very clever at the same time. The trivida® wheel and the innovative triatec technology are subject of the granted European patent EP 3222263.

Technical data

Available in	Plastic PT20	Carbon CC21
Wheel diameter	24"	
Maximum permissible camber	3°	
Load-bearing capacity	Wheelchair 150 kg	Wheelchair 200kg
Axle	Quick-change axle 2 bearings for enhanced stability	
Weight of segment	1050 g / 3.9kg per wheel	700g / 3.0kg per wheel
Design	Spoke wheel, matt surface, black	High-gloss carbon disc, black
Special feature	The robust trivida® wheel in reliable plastic quality offers an optimal price-performance ratio.	The high-strength trivida® carbon wheel in the modern design convinces with its accurate handling and particularly long life.

Guarantee	24 months
Certification	CE ISO13485
Scope of delivery	Quick-change axle: \varnothing 12,7mm OR 12 mm with 50 mm insertion depth, Compatible assembly bush M16 with 50 mm insertion depth bushes for 1/2" axle system Stainless steel nut M16