

Easy GO

Folding Paediatric Wheelchair




Zippie **Easy Go** embodies versatility and efficiency, designed to meet the needs of children during short-term post-surgical recovery or for longer mobility needs. With a simple folding design for ease of transport, the Zippie **Easy Go** is a great companion for families on the go.

Built in Growth

Seat width adjusts from 12-15"

Seat height adjusts from 16-18"

Seat depth adjusts from 14-16"



Parents will appreciate the thoughtful design of easy folding, height adjustable push handles, adjustable armrests and padded adjustable back upholstery



Height adjustable push handles



Centre of gravity adjustable



Foldable push handles with ergonomic design



Adjustable velcro straps



Detachable & height adjustable armrest

Choose Zippie **EasyGo** for unmatched flexibility and comfort!

Specifications: ZEG1214	
Seat width adjustable	12", 13", 14", 15"
Seat height adjustable	16", 17", 18"
Seat depth	14 - 16"
Centre of gravity adjustable	0 - 3"
Stroller handle	Adjustable ergonomic shape
Rear wheels	22" x 1 3/8"
Casters	6" x 1 1/4"
Weight capacity	75 kg
Product weight	19 kg
Warranty	3 Year Frame Warranty