

POWER
MOBILITY

@POWER_MOBILITY

Tuki Catalog

PRODUCT YEAR 2025



Design that moves with you.

At Power Mobility, we're proud to supply Tuki's range of shower commodes—designed and manufactured in Aotearoa New Zealand. With a focus on independence, simplicity, and real-life functionality, Tuki products are built to move with people, not hold them back. We're committed to making assistive technology seamless, so people can focus on living—not the equipment.

@POWER_MOBILITY



Tuki Travel Commode

One of the lightest self-propelled commodes on the market, it packs down simply and intuitively into a wheelie bag that will fit in the car boot, or check it in under your weight allowance for your flight.

Scan QR code below to watch a video



Key Features:

- Designed and made in New Zealand
- High grade anodised aluminium frame (120kg Weight Limit)
- Custom frame dimensions available
- No tools required
- Simple setup & pack-down
- Durable canvas travel wheelie bag
- Flip-over armrests for easy transfers
- Handmade fully padded memory foam seat w/ desired cut-out/bite
- Quick drying canvas backrest
- Side frame cutaways
- Self Propelled wheels and quick-release axles
- Durable and easy to clean design
- Front braking medical grade casters
- Flip-up footplate with adjustable height
- Clearance for positioning over toilet
- Composite scissor lock Brakes



Tuki Tilt Commode

The TUKI Tilt-In-Space Shower Chair offers superior support for the user

Key Features:

- Designed and made in New Zealand
- High grade stainless steel frame
- Swing over armrests for easy transfers
- Handmade fully padded memory foam seat with desired cut-out/bite
- Quick drying canvas backrest
- 0 to 35° posterior tilt guided by braking gas struts
- Adjustable back angle
- Front frame cutaways
- Braked attendant castors
- Durable and easy to clean design
- Drop in swing away foot paddles
- Clearance for positioning over toilet



Tuki Tilt Commode -Collapsible

Now available - The TUKI Collapsible Tilt-In-Space which at a push of a button can collapse into a travel bag.





Tuki Fixed Commode

The TUKI Fixed Frame Shower Chair combines durability and stability with elegance

Key Features:

- Designed and made in New Zealand
- High grade stainless steel frame
- Swing over armrests for easy transfers
- Handmade fully padded memory foam seat with desired cut-out/bite
- Quick drying canvas backrest
- Front and side frame cutaways
- Adjustable backrest angle
- Self propelled wheels and quick-release axles
- Durable and easy to clean design
- Front braking medical grade castors
- Flip-up footplate with adjustable height
- Clearance for positioning over toilet



Tuki Folding Commode

The TUKI Folding Shower Chair combines foldability and stability with elegance.

Key Features:

- Designed and made in New Zealand
- High grade stainless steel frame (140kg weight limit)
- No tools required
- Simple setup & pack-down
- Durable canvas travel wheelie bag
- Swing over armrests for easy transfers
- Handmade fully padded memory foam seat w/ desired cut-out/bite
- Quick drying canvas backrest
- Front and side frame cutaways
- Adjustable backrest angle
- Self Propelled wheels and quick-release axles
- Durable and easy to clean design
- Front braking medical grade casters
- Flip-up footplate with adjustable height
- Clearance for positioning over toilet
- Folded Dimensions (Approx): 470 x 600 x 300mm

POWER
MOBILITY

Contact

@POWER_MOBILITY

Contact

2A/452 Bilsen Rd

Geebung

Queensland

(07) 3265 4663

@power_mobility

www.powermobility.com.au

