

TRAVEL COMMUNE



Key Features:

- High grade anodised aluminum frame
- Flip-over armrests for easy transfers
- Handmade fully padded memory foam seat with desired cut-out/bite
- Quick drying canvas backrest
- Side frame cutaways
- Self propelled wheels and quick-release axles
- Durable and easy to clean design
- Front braking medical grade castors
- Flip-up footplate with adjustable height
- Clearance for positioning over toilet
- Snap on rear attendant castors
- Simple setup & pack-down
- Durable canvas wheelie bag

Specifications:

- * Weight 14kg (In Bag)
- * Weight limit 120kg
- * Overall Dimensions 28.5 x 24.5 x 38" (Standard LxWxH)
- * Custom wheel size (standard 20 or 22")
- * Custom seat width (standard - 18")
- * Custom seat height (standard - 20")
- * Custom seat length (standard - 16")
- * Custom seat rake (standard 1/2")
- * Custom back angle (4 deg standard)

Optional Extras:

- * Backrest: Canvas (standard) or tension adjustable
- * Seat Padding: 1/2 or full gel overlay
- * Footplate: Centre mount flip up (standard), dual swing away or drop in.
- * Brakes: Composite scissor lock.
- * Armrests: Swing over (standard)
- * Light-weight travel bedpan & carrier

Safe use, operation and maintenance.

The Tuki Travel Commode has a maximum user weight of 120kg and is intended for limited use and is NOT designed for everyday use.

Rear wheel tires are solid and require no checking or adjusting of tire pressure

Castor wheel axle/bearings are fully sealed and are maintenance free, lubrication of castor pintles should be carried out by an approved technician.

TURNING CIRCLE:

