

TILT-IN-SPACE COMMmode




Safe use, operation and maintenance.

The Tuki Tilt-In-Space Commmode has a maximum user weight of 140kg and is designed for everyday use.

Castor wheel axle/bearings are fully sealed and are maintenance free, lubrication of castor pintles should be carried out by an approved technician.

TURNING CIRCLE:



Key Features:

- High grade stainless steel frame
- Swing over armrests for easy transfers
- Handmade fully padded memory foam seat with desired cut-out/bite
- Quick drying canvas backrest
- 0 to 35° posterior tilt guided by braking gas struts
- Adjustable back angle
- Front frame cutaways
- Braked attendant castors
- Durable and easy to clean design
- Drop in swing away foot paddles
- Clearance for positioning over toilet

Specifications:

- * Weight 30kg
- * Weight limit 140kg
- * Custom seat width (standard - 18")
- * Custom seat height (standard - 17 1/4")
- * Custom seat length (standard - 16 1/2")
- * Custom back height (standard 17 1/4")
- * Custom back angle (Standard 0 deg)
- * Custom armrest height (Standard 9 1/2")

Optional Extras:

- * Backrest: Canvas (standard), hardshell or tension adjustable
- * Seat Padding: 1/2 or full gel overlay
- * Footplate: Dual swing away
- * Armrests: Swing over (standard), drop-in T style, gutter armrests
- * Headrests & Supports: Thigh Guides, headrest (ball and socket or heavy duty) and swing away/fixed laterals
- * Bedpan & carrier