

QUANTUM® **EDGE® 4**

A QUANTUM LEAP.



FEATURING 9 NEW
COLOR OPTIONS!



THE MOST ADVANCED EDGE YET

The Edge® 4 is our most advanced model yet. With 180 mm (7") front and rear casters and powerful 4-pole motors, it delivers speeds up to 10 km/h and confidently tackles everyday obstacles. A new battery tray allows easier access and supports 75 Ah batteries for extended range. LED accent lights are integrated into the casters and fenders for greater visibility and style.

Enhanced suspension and damping smooth out rough terrain, giving you comfort, control, and performance wherever you go..



QUANTUM®



quantumrehab.com.au

Specifications

DRIVE WHEELS⁹	35.5cm (14") pneumatic or flat free tyres Optional: 35.5cm (14") Community flat free tyres available
CASTER WHEELS	17.8cm flat free (7")
MAXIMUM SPEED¹	Up to 10 kph (7.7 km/h with iLevel®)
RANGE PER CHARGE^{1,8}	Up to 29.1 km
GROUND CLEARANCE	7.28cm (2.87")
TURNING RADIUS	57.9cm (22.8")
LENGTH	97.2cm (38.27")
BASE WIDTH	61cm (24.02")
SEATING SIZES	
TRU Balance® 4 Anterior Tilt⁵	Width: 40-55.8cm (16-22") Depth: 35.5-55.8cm (16-22")
TRU Balance® 3 Static Seat	Width: 30.5cm-55.8cm (12-22") Depth: 30.5cm-50.8cm (12-20")
TRU Balance® 3 Power Tilt	Width: 30.5cm-55.8cm (12-22") Depth: 30.5cm-55.8cm (12-22")
TRU Balance® 3 Tilt and iLevel®	Width: 30.5cm-55.8cm (12-22") Depth: 30.5cm-55.8cm (12-22")
SEAT-TO-FLOOR HEIGHTS	
TRU Balance® 4 Anterior Tilt⁵	46.3cm-53.9cm (18.25-21.25")
TRU Balance® 3 Power Tilt⁵	45cm-50.1cm (17.75-19.75")
TRU Balance® 3 Tilt & iLevel®⁵	44.45cm-49.5cm (17.5-19.5")
BATTERY SIZE³ MK 70Ah(2)	
BATTERY WEIGHT⁷ 23.3kg each	
AVAILABLE ELECTRONICS⁶ 120A Q-Logic 3	
BATTERY CHARGER 8A Off-board	
MOTOR PACKAGES 4-Pole	
WEIGHT CAPACITY⁴ 150 kg.	
BASE WEIGHT 193 kg	

1) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 150 kg. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) All Quantum Rehab® units can be ordered with specialty controls.
 7) Battery weight may vary +/- 1.5kg. based on manufacturer.
 8) Range shown is measured using Group 24 batteries
 9) Pneumatic not available in community tyre option.

Features

- Independent SRS (Smooth Ride Suspension)
- Standard led front castor & rear fender lighting package
- USB charger included
- Wide array of color choices
- Q-Logic 3 Drive Controls
- Made in Australia of foreign and domestic parts

Options

- TRU-Balance® 4 Power positioning: 10 degree or 30 degree anterior tilt Options with memory seating
- TRU-Balance® 3 Power Positioning: iLevel® option allows 30.5 cm (12") of seat elevation while at walking speed (7.6 km/h) Power Tilt, Recline, and Elevating Legrests that can power to floor for easy transfers.
- Accu-Trac advanced tracking technology
- ISO 7176-19 Occupied Crash-Test Approved
- NEW backup camera
- Additional options available on scrips / configurator

Colors

